**RAPID HTA REVIEW**

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GORE® Cardioform ASD Occluder experience in transcatheter closure of "complex" atrial septal defects.

Santoro G, Cuman M, Pizzuto A, Haxhiademi D, Lunardini A, Franchi E, Marrone C, Pak V, Assanta N, Cantinotti M.

OBJECTIVE: To evaluate the GORE® Cardioform ASD Occluder (GCO) (WL Gore & Associates, Flagstaff, AZ) device for "complex" atrial septal defects (ASD) closure.

BACKGROUND: Transcatheter ASD closure is still challenging in "complex" clinical/anatomic settings. This study evaluated the results of GCO in closure of "complex" ASD in a tertiary referral center.

METHODS: Between January 2020 and March 2021, 72 patients with significant ASD were submitted to transcatheter closure with GCO at our Institution. Based on clinical/anatomic characteristics, they were classified as "complex" (n = 36, Group I) or "simple" (n = 36, Group II). We considered as "complex", defects with rim deficiency (< 5 mm) other than antero-superior, relatively large (diameter/patient weight > 1.2 or diameter/patient BSA > 20 mm/m²) or within a multifenestrated septum. Procedure results and early outcome were compared

between the groups.

RESULTS: Absolute and relative ASD size (20 ± 4 vs. 15 ± 3 mm, $p < 0.0001$; 0.9 ± 0.3 vs. 0.4 ± 0.2 mm/kg, $p < 0.0001$; 23 ± 7 vs. 12 ± 5 mm/m², $p < 0.0001$), QP/QS (2.0 ± 0.8 vs. 1.4 ± 0.3 , $p < 0.001$), procedure and fluoroscopy times (73 ± 36 vs. 43 ± 21 min, $p < 0.0001$; 16 ± 9 vs. 9 ± 4 min, $p < 0.0001$), procedure feasibility (94.4 vs. 100%, $p < 0.0001$) and overall complication rate (13.9 vs. 0%, $p < 0.0001$) were significantly different between the groups. Successful closure of "surgical" ASDs was achieved in 92% of cases. Complete closure at last follow-up evaluation did not significantly differ between the groups (97.1 vs. 100%, $p = \text{NS}$), as was wireframe fractures rate (49.1% in the overall population), without clinical, EKG and echocardiographic consequences.

CONCLUSIONS: Percutaneous treatment with GCO device is effective and safe in high percentage of "complex" ASDs.

2. Int J Cardiol. 2021 Mar 15;327:68-73. doi: 10.1016/j.ijcard.2020.11.029. Epub 2020 Nov 18.

Trans-catheter atrial septal defect closure with the new GORE® Cardioform ASD occluder: First European experience.

Santoro G, Castaldi B, Cuman M, Di Candia A, Pizzuto A, Sirico D, Cantinotti M, Garibaldi S, Pak V, Di Salvo G.

BACKGROUND: This perspective, observational study evaluated safety and efficacy of the GORE® Cardioform ASD Occluder (WL Gore & Associates, Flagstaff, AZ), compliant and potentially innovative prosthesis recently approved for closure of ostium secundum atrial septal defects (ASD).

METHODS: Between January and June 2020, 43 unselected patients with -significant ASD were submitted to trans-catheter closure with GORE® Cardioform ASD Occluder at two high-volume Italian Pediatric Cardiology centers. Primary endpoints were procedural success and safety. Secondary endpoints were closure rate and clinical safety at 1-month follow-up.

RESULTS: Patients' age and weight were 8.2 ± 3.9 years (range 3-21, median 9.9) and 29.6 ± 15.3 kg (range 16-57, median 33.3), respectively. ASD diameter was 16.6 ± 4.5 mm (median 10), resulting in QP/QS of 1.7 ± 0.7 (median 1.6). Seventeen pts. (39.5%) were considered "surgical" candidates due to challenging septum morphology, ASD rim deficiency or ASD diameter/patient weight ratio ≥ 1.2 . Device placement was successfully achieved in all but one patient (97.7%), in whom it embolized early after deployment, resulting in rescue surgical repair. No cross-over with different devices was recorded. Median procedure and fluoroscopy times were 40 and 6.8 min, respectively. Major adverse events were recorded in 7.0% (3 pts). Complete closure rate was 78.5% at

discharge, rising to 92.9% (39/42 pts) at 1 month evaluation, without cardiac or extra-cardiac adverse events. "Challenging" procedures were more time-consuming but as effective and safe as the "simple" ones.

CONCLUSIONS: The GORE® Cardioform ASD Occluder device was highly effective and versatile in closure of ASDs with different anatomy and size, even in challenging settings.

3. Catheter Cardiovasc Interv. 2017 Sep 1;90(3):495-503. doi: 10.1002/ccd.26907. Epub 2017 Jan 27.

Initial clinical experience with the GORE®) CARDIOFORM ASD occluder for transcatheter atrial septal defect closure.

de Hemptinne Q, Horlick EM, Osten MD, Millán X, Tadros VX, Pighi M, Gonzalez Barlatey F, Alnasser SM, Miró J, Asgar AW, Ibrahim R.

OBJECTIVES: To report the initial clinical experience with a novel atrial septal defect (ASD) closure device, the GORE® CARDIOFORM ASD Occluder (GCO).

BACKGROUND: Transcatheter closure has become the treatment of choice for secundum ASD. A wide range of occluder devices are available, but concern has been raised about the risk of cardiac erosion associated with rigid devices and future access to the left atrium. **METHODS:** Retrospective chart review of patients treated with the GCO at three Canadian centers. Primary outcomes were procedural success and residual shunting at follow-up, as well as 30-day major adverse events. Secondary outcomes included new onset atrial arrhythmias, wire frame fractures (WFF), and all cause mortality. Clinical, echocardiographic, procedural data, and follow-up outcome variables were collected in each participating hospital.

RESULTS: Between February and December 2015, 26 patients (5 children and 21 adults) underwent transcatheter ASD closure with the GCO and were included in the study cohort. Procedural success was achieved in 22 of 26 patients (85%) and no major procedural complications were observed. Two patients (8%) presented new

onset atrial tachyarrhythmia during early follow-up (0-30 days). Follow-up echocardiography (median of 119 days [IQR: 92-146]) demonstrated no residual shunt in all implanted patients. After a median clinical follow-up of 174 days (IQR: 135-239), one patient died of an unrelated cause, there were no documented major adverse cardiovascular events. Fluoroscopic imaging of the device was performed in 20 patients (91%), and WFF was noted in five cases. No clinical consequence or device dysfunction was observed in these patients.

CONCLUSIONS: In this first-in-man multicenter study, the GCO was safe and effective for ASD closure, with no major adverse events or residual shunt at follow-up. © 2017 Wiley Periodicals, Inc.

4. Catheter Cardiovasc Interv. 2020 Jun 1;95(7):1285-1295. doi: 10.1002/ccd.28728. Epub 2020 Jan 14.

ASSURED clinical study: New GORE® CARDIOFORM ASD occluder for transcatheter closure of atrial septal defect.

Sommer RJ, Love BA, Paolillo JA, Gray RG, Goldstein BH, Morgan GJ, Gillespie MJ; ASSURED Investigators.

OBJECTIVE: Assess clinical performance of a new device for transcatheter closure of atrial septal defect (ASD).

BACKGROUND: Previously-approved ASD Closure devices have known limitations. Device erosion has been associated with the AMPLATZER® septal occluder in patients with retro-aortic rim deficiency (<5 mm), while defects ≥18 mm are too large for the GORE® CARDIOFORM septal occluder. The GORE® CARDIOFORM ASD occluder (GCA), a hybrid of the approved devices, was designed to expand the eligible ASD population.

METHODS: One-hundred and twenty-five ASD patients were enrolled in a prospective, multicenter registry. Descriptive clinical endpoints included: technical implant success, closure success, serious adverse events (SAE), clinically significant new arrhythmia, and wire frame fracture. Procedural outcomes and adverse events were adjudicated by an Echocardiography Core Lab and Independent Data Review Board, respectively.

RESULTS: Median subject age was 12.3 years (range 2.9-84.7), with 72% of patients ≤18 years old. Median ASD stop-flow diameter was 17.0 mm (8.0-30.0), with 43% ≥18 mm. Deficient retro-aortic rim occurred in 57% of subjects, and 30% had both diameter ≥ 18 mm and deficient rim. Technical Implant Success was achieved in 120/125 (96%), though three devices were removed within 24 hr. At the scheduled 6-month evaluation, 112/117 returned for evaluation. All 112 had closure success. SAE occurred in 6/125 (4.8%) and 6/125 (4.8%) had clinically significant new arrhythmia. Wire frame fractures occurred in 37/104 (36%), without any associated clinical sequelae, residual shunt, or device instability.

CONCLUSION: The GCA performed effectively and safely in this initial cohort, and led to FDA approval based on this data.

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