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Alice Corzani

**BIBLIOGRAFIA**

1. Axiostat, Instructions for use, indirizzo internet, <https://axiobio.com/axiostat-haemostatic-dressings/>
2. Minici R, Serra R, Maglia C, Guzzardi G, Spinetta M, Fontana F, Venturini M, Laganà D. Efficacy and Safety of Axiostat® Hemostatic Dressing in Aiding Manual Compression Closure of the Femoral Arterial Access Site in Patients Undergoing Endovascular Treatments: A Preliminary Clinical Experience in Two Centers. J Pers Med. 2023 May 11;13(5):812. doi: 10.3390/jpm13050812. PMID: 37240982; PMCID: PMC10220903.
3. Nakache A, Darmon A, Molho A, Steinecker M, Nejjari M, Digne F. Evaluation of the safety and efficacy of the Axiostat® dressing device to achieve radial artery access hemostasis: The R3A study. Catheter Cardiovasc Interv. 2024 Jun 16. doi: 10.1002/ccd.31129. Epub ahead of print. PMID: 38881025.

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RAPID HTA REPORT		
No. of the request	Date of request	Applicant
323	10/06/2024	Diagnostic and therapeutic neuroradiology - AOUS
Type of report		
New report		X
Update of a previous report		
If update, state the reason:		

General technology data			
Trade name			
AXIOSTAT VASCULAR			
Generic name			
Chitosan dressing to control bleeding from vascular access sites			
Manufacturer name			
AXIO BIOSOLUTIONS PVT.LTD			
Supplier name			
SVAS BIOSANA S.p.A.,			
RDM	REF		
1704877	V55 (5x5x1 cm), V35 (3.5x3.5x1 cm), V25 (2.5x2.5x1 cm)		
Type	CE mark (date)	Risk class	FDA approval
1	03/03/2020 (first 02/07/2018)	III	Yes
CND			
M040599 - Hemostatic dressings - other			
Clinical problem and rationale for the request			
Pediatric postprocedural hemostasis in femoral punctures as young patients are at risk of introducer insertion site bleeding in the early moments and within 12 hours postprocedure.			
Indications for use			
<p>The data sheet directions are given below [1]:</p> <ul style="list-style-type: none"> <li>-VASCULAR USE Axiostat Vascular Hemostatic Dressing is a sterile, single-use, nonabsorbable hemostatic dressing intended to control bleeding from vascular access sites by reducing ambulatory time in patients undergoing diagnostic or interventional endovascular procedures using cannulated guidewire surgical procedures for peripheral access.</li> <li>-USE IN EMERGENCY Axiostat is a sterile, single-use, nonabsorbent hemostatic dressing intended for use in the temporary control of bleeding wounds. Axiostat blocks moderate to severe bleeding due to cuts, abrasions, lacerations, venous or arterial bleeding.</li> <li>-MILITARY USE Axiostat is a sterile, single-use, nonabsorbent hemostatic dressing intended for use in the temporary control of bleeding wounds. Axiostat blocks moderate to severe bleeding from cuts, abrasions, lacerations, venous</li> </ul>			

**General technology data**

or arterial bleeding. -DENTAL USE Axiostat dressing is a sterile, nonabsorbable hemostatic dressing intended to control oral bleeding during dental procedures.
<b>Target patient</b>
Pediatric patients following procedures that involved femoral puncture.
<b>Main competitors</b>
Fibrillar Tabotamp (now Surgicel).

**Technology details**

<b>Description</b>
According to the data sheet, the Axiostat device is a hemostatic dressing composed of 100% pure, sterile, single-use, nonabsorbable chitosan based biomaterial for the control of bleeding wounds. The hemostatic agent acts through an interaction between the erythrocyte cell membrane (negative charge) and the protonated biomaterial (positive charge) leading to platelet involvement and rapid thrombus formation. The biomaterial, chitosan used for this device, is derived from crustaceans. The device is packaged in moisture-proof containers and terminally sterilized by gamma radiation.
<b>Elements of innovation</b>
The clinician reports that compared with other hemostatic dressings awarded to tender, AXIOSTAT provides the following advantages: improved water and blood absorption and adhesion capacity, faster hemostasis due to the intrinsic cationic charge of the re-engineered chitosan. As there is no literature present these innovative/improving features are assumed and potential.

**Clinical and economic evidence**

<b>Clinical studies</b>
The Pubmed database search with keyword AXIOSTAT conducted on 06/25/2024 yielded 8 results, of which only one was relevant to the subject of the analysis [2]. The objective of that study (Minici et al. [2]) was to evaluate the safety and efficacy of Axiostat in facilitating closure by manual compression of the femoral arterial access site in patients undergoing endovascular treatment. Specifically, the study retrospectively evaluated 120 adult patients (age under 18 years was reason for exclusion) who had undergone Axiostat-assisted manual compression of the femoral arterial access site after endovascular treatments from July 2022 to February 2023. The primary efficacy endpoint was the primary technical success rate achieved with manual compression assisted by Axiostat hemostatic dressing. Primary technical success was defined by adequate hemostasis achieved with 7 minutes (in cases of 4-6 Fr size introducers) or 10 minutes (in cases of 7-8 Fr size introducers) of manual compression aided by the Axiostat hemostatic dressing. Secondary technical success was defined as adequate hemostasis after the above scheme or after an additional period of Axiostat-assisted manual compression. The primary safety endpoint was clinical success rate. Clinical success was defined by the absence of bleeding-related VASC (vascular access site complications) at follow-up. Secondary endpoints included secondary technical success rate, TTH (time required to achieve hemostasis), TTA (time after which ambulation is possible), and rate of VASCs (complications at the vascular access site). The secondary objective of the study was to evaluate and compare efficacy and safety between two subgroups (Group 1: introducer size 4-6 Fr; Group 2: introducer size 7-8 Fr). A total of 120 patients were identified. Primary technical success was achieved in 110 (91.7%) patients, with adequate hemostasis achieved in all cases of prolonged manual compression requirement. The mean TTH and TTA were 8.9 (±3.9) and 462 (±199) minutes, respectively. Clinical success was achieved in 113 (94.2%) patients, with bleeding-related VASCs observed in 7 (5.8%) patients.

The rate of vascular access site complications (VASCs) was 8.3%, related to two (1.7%) cases of hematoma, five (4.1%) pseudoaneurysms and three (2.5%) cases of dissection. No major hemorrhagic events occurred. According to the CIRSE complication classification system, only grade 2 or 3 events occurred. According to SVS reporting standards, no major VASCs (grade 3) were detected. Only one case of pseudoaneurysm required an invasive procedure (percutaneous thrombin injection). Medical treatment consisted of prolonged application of a pressure dressing in two cases (1.7%), ultrasound-guided compression of pseudoaneurysms in four cases (3.3%), and administration of antiplatelet therapy to prevent complications of dissections in three cases (2.5%).

A comparison of data was performed between patients undergoing endovascular treatment with introducer sheaths between 4 and 6 Fr in size (Group 1) and patients undergoing endovascular treatment with introducer sheaths between 7 and 8 Fr in size (Group 2).

No statistically significant differences were observed between the group undergoing endovascular treatment with introducer sheaths between 4 and 6 Fr and the group undergoing endovascular treatment with introducer sheaths between 7 and 8 Fr in terms of age (69.2 vs. 65.5 years;  $p = 0.317$ ), coagulopathy (11 vs. 1 patient;  $p = 0.291$ ), INR (1.2 vs. 1.2;  $p = 0.937$ ), antiplatelet or anticoagulant therapy (51 vs. 19 patients;  $p = 0.343$ ), systolic blood pressure  $\geq 180$  mmHg (10 vs. 6 patients;  $p = 0.201$ ), VASC (9 vs. 1;  $p = 0.450$ ) and bleeding-related VASC (6 vs. 1;  $p = 1$ ). The number of patients who received an intraoperative thrombolytic agent was significantly higher in the 7-8 Fr group (25% vs. 1.1%;  $p < 0.001$ ). The primary technical success rate was significantly higher in the 4-6 Fr group (97.8% vs. 71.4%;  $p < 0.001$ ).

The study concludes that the results of the retrospective investigation demonstrate that manual compression aided by Axiostat hemostatic dressing is effective and safe for achieving femoral arterial access site hemostasis in patients undergoing endovascular treatment with a 4-8 Fr introducer sheath. The absence of a control group limits the strength of the investigation. Larger, prospective, randomized, controlled studies are desirable to confirm these preliminary data and compare the efficacy and safety profiles of Axiostat-assisted manual compression, manual compression alone, and VCDs (vascular closure devices).

In addition, coagulopathy and systolic blood pressure recorded before sheath removal are reliable and independent predictors of the occurrence of bleeding-related VASCs.

Also of note is the study by Nakache et al. [3], although not relevant by indication for use because it related to achieving hemostasis in patients undergoing transradial coronary angioplasty, because according to it Axiostat had been shown to be effective in achieving hemostasis after TRCA, but its use was found to be associated with an unexpectedly high incidence of radial thrombosis (in fact, the study was prematurely terminated after the inclusion of 41 patients instead of the expected 150).

Pubmed search (<https://pubmed.ncbi.nlm.nih.gov/>): keyword: axiostat; search date: 06/25/2024.

**Ongoing clinical trials**

No studies found on clinicaltrials.gov

**Guidelines**

No guidelines

**Cost-effectiveness analysis**

None

**HTA Report**

None

**Expected benefits**

The reported advantages should be: improved water and blood absorption capacity, faster hemostasis due to the intrinsic cationic charge of the reengineered chitosan.

As there is no literature related to this DM these innovative/improving features are assumed and potential.

Price and cost of therapy per patient			
Product (Manufacturer)	Unit price (euro)	Cost of therapy per patient (euro)	Annual requirement (No. pieces)
Axiostat vascular 5x5x1 cm (AXIO BIOSOLUTIONS)	80 per envelope	80 per envelope	120

Price and cost of therapy per patient with treatment alternatives already in use		
Product (Manufacturer)	Unit price (euro)	Cost of therapy per patient (euro)
Surgicel 2.5x5.1 cm (Ethicon)	31.20 per envelope	31.20 per envelope

Economic and organizational impact
Not evaluated due to lack of specific information.

Reimbursability information			
Principal diagnosis (ICD9-CM code).	Intervention (ICD9-CM Code )	DRG Code	Rate (euro)
19.05	99.25	variable	-

Evaluation of innovativeness (according to Regional Resolution No. 737/2022, <a href="#">link</a> and No. 1244/2022 <a href="#">link</a> )	
Innovative device (Y/N)	N
If yes, please indicate which Criteria 1, 2 and 3 are met:	

Summary data		
Request number	Date of request	Applicant
323	10/06/2023	Diagnostic and therapeutic neuroradiology - AOUS
Technology		
Axiostat Vascular, chitosan dressing intended to control bleeding from vascular access sites by reducing ambulatory time in patients undergoing diagnostic or interventional endovascular procedures using cannulated guidewire surgical procedures for peripheral access.		
Conclusions		
Axiostat is a chitosan dressing whose efficacy and safety profile remains to be defined since to date this dressing has only been evaluated in a single retrospective single-arm trial conducted in 120 adult patients (age under 18 years was a reason for exclusion from the study) [2]. To define the real clinical benefit, comparative studies toward manual compression alone or toward manual compression assisted with other dressing (Surgicel) are needed. Axiostat also has a higher price than Surgicel. The opinion is therefore unfavorable.		
Date of report writing		
20/09/2024		
Author(s) of the report		
Beatrice Tognarelli		
Referring pharmacist		
Alice Corzani		

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*Centro operativo*

## BIBLIOGRAPHY

1. Axiostat, Instructions for use, Internet address, <https://axiobio.com/axiostat-haemostatic-dressings/>
2. Minici R, Serra R, Maglia C, Guzzardi G, Spinetta M, Fontana F, Venturini M, Laganà D. Efficacy and Safety of Axiostat® Hemostatic Dressing in Aiding Manual Compression Closure of the Femoral Arterial Access Site in Patients Undergoing Endovascular Treatments: A Preliminary Clinical Experience in Two Centers. J Pers Med. 2023 May 11;13(5):812. doi: 10.3390/jpm13050812. PMID: 37240982; PMCID: PMC10220903.
3. Nakache A, Darmon A, Molho A, Steinecker M, Nejjari M, Digne F. Evaluation of the safety and efficacy of the Axiostat® dressing device to achieve radial artery access hemostasis: The R3A study. Catheter Cardiovasc Interv. 2024 Jun 16. doi: 10.1002/ccd.31129. Epub ahead of print. PMID: 38881025.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Drafting of the report by the Operations Center, Regional Decree No.17610 of September 7, 2022.

For more information email [centro.operativo.htart@regione.toscana.it](mailto:centro.operativo.htart@regione.toscana.it)