
**RAPID HTA REVIEW**

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Caterina Sgromo

BIBLIOGRAFIA

1. Zhang Y, Xiang D, Lu Q, Wu M, Cui J. A systematic review and meta-analysis of the performance of flowdiverting stents in the treatment of peripheral and visceral artery aneurysms. *Catheter Cardiovasc Interv.* 2021 Feb 15;97(3):461-469. doi: 10.1002/ccd.29373. Epub 2020 Nov 11. PMID: 33175422.
2. Rabuffi P, Bruni A, Antonuccio EGM, Ambrogi C, Vagnarelli S. Treatment of visceral artery aneurysms and pseudoaneurysms with the use of cerebral flow diverting stents: initial experience. *CVIR Endovasc.* 2020 Sep 13;3(1):48. doi: 10.1186/s42155-020-00137-y. PMID: 32886269; PMCID: PMC7474014.
3. Giorgakis E, Chong B, Oklu R, Jaroszewski DE, Knuttinen G, Mathur AK. Successful treatment of visceral pseudoaneurysm after pancreatectomy using flow-diverting stent device. *Ann Hepatobiliary Pancreat Surg.* 2020 Feb;24(1):114-118. doi: 10.14701/ahbps.2020.24.1.114. Epub 2020 Feb 27. PMID: 32181440; PMCID: PMC7061036.
4. Eldem G, Erdoğan E, Peynircioğlu B, Arat A, Balkancı F. Endovascular treatment of true renal artery aneurysms: a single center experience. *Diagn Interv Radiol.* 2019 Jan;25(1):62-70. doi: 10.5152/dir.2018.17354. PMID: 30272561; PMCID: PMC6339632.
5. Jamshidi M, Rajabian M, Avery MB, Sundararaj U, Ronsky J, Belanger B, Wong JH, Mitha AP. A novel selfexpanding primarily bioabsorbable braided flow-diverting stent for aneurysms: initial safety results. *J Neurointerv Surg.* 2020 Jul;12(7):700-705. doi: 10.1136/neurintsurg-2019-015555. Epub 2019 Nov 27. PMID: 31776171.

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
RAPID HTA REVIEW			
Request No.	Date of request	Applicant	
285	05/04/2023	AOUC - Interventional Radiology	
Type of report			
New report			S
Updating a previous report			No
If updating, please state the reason:			
General technology data			
Trade name			
DERIVO PERIPHERAL EMBOLISATION DEVICE			
Generic name			
Peripheral flow-diversion stent			
Manufacturer name			
Acandis GmbH			
Supplier name			
CrossMed			
RDM	REF		
2290762; 290070; 2290772; 2290773; 290777; 2290779; 2290780; 2290781; 2290782; 2290786; 2290787; 2290788; 2290789; 2290790; 2290791; 2290795; 2290796; 2290797; 2290798; 2290799; 2290703; 2290804; 2290805; 2290806; 2290807; 2290811; 2290812.	da 01-102310 a 01102313; 01-102363; 01-102364; da 01-102410 a 01-102413; 01-102463; 01-102464; da 01-102511 a 01102513; 01- 102563; 01102564; da 01-102514 a 01-102516; 01- 102568; 01- 102569; 01-102520; 0-102521; 01-102573; 01102574.		
Type	CE marking (date)	Risk class	FDA approval
1	516802 MR2 (05/2021)	IIb	No
CND			
P07040202 (COATED PERIPHERAL VASCULAR STENTS)			
Clinical problem and rationale for the request			
Treatment of visceral and peripheral aneurysms.			

General technology data
Indication of use
Derivo Peripheral Embolisation Device is a flow-diversion stent indicated for the treatment of peripheral and visceral arterial aneurysms (visceral arteries: spleen, liver, inferior/superior mesenteric artery, celiac, renal; peripheral arteries: popliteal, subclavian, iliac and femoral arteries) that cannot be treated with other techniques endovascular or where other endovascular or surgical techniques present a higher risk of treatment.
Target patient
Patients with visceral and peripheral aneurysms that cannot be treated with other endovascular techniques or where other endovascular or surgical techniques present a higher risk of treatment.
Main competitors
Spirals and embolising fluid used additionally. The flow-diverting stents in use today (e.g. Fred, Silk, Derivo) are indicated for the treatment of intracranial aneurysms and are therefore not approved for use on visceral vessels or peripheral arteries.
Target patient
Patients with visceral and peripheral aneurysms that cannot be treated with other endovascular techniques or where other endovascular or surgical techniques present a higher risk of treatment.
Technological details
Description
Derivo peripheral Embolisation Device is a self-expanding device, braided with nitinol (nickel-titanium alloy) wires with a platinum core, pre-loaded on a transport wire in an introducer. Derivo peripheral Embolisation Device is designed to be introduced via a microcatheter with an internal diameter of 0.027" or 0.039", a depending on size.
Elements of innovation
Derivo Peripheral Embolisation Device is the only flow-diversion stent with an indication for use on visceral vessels or peripheral arteries. Its innovative features are its small size, degree of flexibility and the possibility to be loaded on microcatheters, which allow the Derivo stent to be routed further downstream reducing the risk of ischaemia. Flow diversion allows the collateral branches (efferent vessels) to be kept open while excluding the aneurysmal sac.
Clinical and economic evidence
Clinical studies
The literature search did not reveal any studies that evaluated the use of the Derivo peripheral embolisation stent for the treatment of peripheral and visceral arterial aneurysms. All studies on the Derivo stent concern the treatment of intracranial aneurysms. Other flow-diversion stents have been studied for peripheral and visceral use, however, the studies are few and mostly case reports. Such stents have been used off-label since their approval is for the treatment of cerebral aneurysms only [1-5].
MEDLINE search PubMed version (site: www.pubmed.org) conducted on 2 May 2023 with the keyword 'Derivo stent'.
Ongoing clinical trials
On ClinicalTrials.gov (https://clinicaltrials.gov/ , keyword "derivative foripher", date 02/05/2023) is not registered no study.
Guidelines

Price and cost of therapy per patient			
Product (Manufacturer)	Price (euro)	Therapy cost per patient (euro)	Annual requirement (N, pieces)
Derivo Embolization device	8,400 (4% VAT excluded)	8.400	18
Cost-effectiveness analysis			
None			
HTA Report			
Price and cost of therapy per patient with existing treatment alternatives			
Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)	
Embolising spirals	900	5.400*	
Onyx embolising fluid	890	890	
		6.290	
*obtained by considering the use of 6 spirals.			
Economic and organisational impact			
The use of the Derivo peripheral stent entails an additional cost of approximately €2,000 per patient.			
Reimbursement procedure related to the use of the requested medical device			
Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM code)	DRG code	Tariff (euro)
747.6 (Other abnormalities of the peripheral vascular system)	39.79 (Repair of other vessels)	555 (Percutaneous cardiovascular interventions with a major cardiovascular diagnosis)	6.393
Assessment of innovativeness (according to Regional Resolution N° 737/2022, link and N° 1244/2022 link)			
Innovative device (Y/N)		N	
If yes, please indicate which Criteria 1, 2 and 3 are met:			
Summary data			
Request number	Date requested	Applicant	
285	5 April 2023	AOUC – Radiologia interventistica	

Summary data
Technology under evaluation
Peripheral Embolisation Device (peripheral flow diversion stent)
Any external experts involved
-
Conclusions
Given that Derivo peripheral embolisation device is the only flow-diversion stent on the market with an indication on visceral vessels or peripheral arteries (all other flow-diversion stents are in fact approved for use in the brain), the CO expresses a favourable opinion on the purchase. However, since, to date, no studies are available to support the efficacy and safety of this device, the CO suggests that the requesting clinician undertake an observational clinical study. The CO is willing, if necessary, to collaborate with the clinician to finalise all the material necessary to submit the study protocol to the ethics committee.
Date of drafting
03/05/2023
Author
Sabrina Trippoli
Pharmacist in charge of request
Caterina Sgromo

BIBLIOGRAPHY

1. Zhang Y, Xiang D, Lu Q, Wu M, Cui J. A systematic review and meta-analysis of the performance of flow-diverting stents in the treatment of peripheral and visceral artery aneurysms. *Catheter Cardiovasc Interv*. 2021 Feb 15;97(3):461-469. doi: 10.1002/ccd.29373. Epub 2020 Nov 11. PMID: 33175422.
2. Rabuffi P, Bruni A, Antonuccio EGM, Ambrogi C, Vagnarelli S. Treatment of visceral artery aneurysms and pseudoaneurysms with the use of cerebral flow diverting stents: initial experience. *CVIR Endovasc*. 2020 Sep 13;3(1):48. doi: 10.1186/s42155-020-00137-y. PMID: 32886269; PMCID: PMC7474014.
3. Giorgakis E, Chong B, Oklu R, Jaroszewski DE, Knuttinen G, Mathur AK. Successful treatment of visceral pseudoaneurysm after pancreatectomy using flow-diverting stent device. *Ann Hepatobiliary Pancreat Surg*. 2020 Feb;24(1):114-118. doi: 10.14701/ahbps.2020.24.1.114. Epub 2020 Feb 27. PMID: 32181440; PMCID: PMC7061036.
4. Eldem G, Erdoğan E, Peynircioğlu B, Arat A, Balkancı F. Endovascular treatment of true renal artery aneurysms: a single center experience. *Diagn Interv Radiol*. 2019 Jan;25(1):62-70. doi: 10.5152/dir.2018.17354. PMID: 30272561; PMCID: PMC6339632.
5. Jamshidi M, Rajabian M, Avery MB, Sundararaj U, Ronsky J, Belanger B, Wong JH, Mitha AP. A novel selfexpanding primarily bioabsorbable braided flow-diverting stent for aneurysms: initial safety results. *J Neurointerv Surg*. 2020 Jul;12(7):700-705. doi: 10.1136/neurintsurg-2019-015555. Epub 2019 Nov 27. PMID: 31776171.

A copy of the document can be downloaded from the website <http://www.regione.toscana.it/-/prodotti-hta>.

Edited by the Operations Centre, established by Decree No. 17610 of 7 September 2022.

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Centro operativo

For further information please write an email to: centro.operativo.htart@regione.toscana.it