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RAPID HTA REPORT		
Request No.	Date of request	Applicant
286	05/04/2023	Urology - AOUC
Type of report		
New report		Yes
Updating a previous report		
If updating, please state the reason:		

General technology data				
Trade name				
ZEPHYR ZSI 475 THREE-COMPONENT PENILE PROSTHESIS				
Generic name				
INFLATABLE FTM THREE-COMPONENT PENILE PROSTHESIS				
Manufacturer name				
ZEPHYR SURGICAL IMPLANTS				
Supplier name				
ZEPHYR SURGICAL IMPLANTS				
RDM	REF			
642556	ZSI 475			
Type	CE marking (date)	Risk class	FDA approval	
1	N° 15749 (05/03/2021)	III	No	
CND				
P080202 - PENILE PROSTHESES				
Clinical problem and rationale for the request				
<p>Phalloplasty is the penis construction operation for women who want to change sex. Phalloplasty addresses two needs, which are to urinate in an erect position and to shower in promiscuous environments. About 50 per cent of women who undergo the phalloplasty procedure then decide to proceed with the implantation of a penile prosthesis for erection.</p> <p>Until recently, all penile prostheses were manufactured for erectile dysfunction in biological males. A few years ago, however, an erectile device was designed specifically for female-to-male transsexuals: the Zephyr ZSI 475 FTM. This prosthesis has a realistic gland shape, a single cylinder, and a specific titanium and silicone anchoring device to anchor the system to the patient's pubic bone. These specific features could solve the problems related to implanting erectile devices in the neofallus and hopefully lead to better surgical results in these patients.</p>				
Indication of use				
<p>The ZSI 475 FTM inflatable penile prosthesis is a device that is inserted into the phalloplasty and allows patients to achieve an erection. This hydraulic implant consists of an inflatable cylinder inserted into the penis, connected to a reservoir containing saline implanted in the pelvis and a hand pump implanted in the scrotum. The prosthesis is mainly made of silicone. The prosthesis is filled under pressure to produce penile rigidity.</p>				
Target patient				
FtM (Female to Male) transgender patients undergoing RCS (surgical sex reassignment).				

General technology data

Main competitors

Penile prostheses (AMS 700) intended for erectile dysfunction that cannot be used in transgender people are currently approved under the same CND.

Technological details

Description

The ZSI 475 FTM Inflatable Penile Implant is a device that is inserted inside the phalloplasty and allows erection. The FTM Inflatable Penile Implant is a device consisting of an inflatable cylinder, connected to a reservoir filled with saline solution, and a manual subcutaneous pump. Four lengths are available: 120, 150, 180 and 210 mm. The inflatable cylinder is inserted into the phalloplasty through a lateral incision. The pump is placed in the scrotum and the reservoir in the pelvis or abdomen. Inflating the implant gives rigidity to the phalloplasty.

Elements of Innovation

The ZSI 475 FTM phalloplasty penile prostheses from the company Zephyr Surgical Implants are, to date, the only penile prostheses on the market specifically designed and developed for gynoandroid surgery.

Clinical and economic evidence

Clinical studies

ZSI 475 FTM phalloplasty penile prostheses were evaluated in two clinical studies [1-2].

The first study was published by Neuville et al. in 2019 [1]. That study included 20 patients with gender dysphoria who underwent surgery for female-to-male transformation. The average age was 37.9 years. Complications after 21 implants included 2 (9.5%) medically treated infections, 1 (4.7%) infection treated with explantation, 2 (9.5%) mechanical failure and 1 (4.7%) malposition. The average follow-up was 8.9 months. 12 patients (85.7%) had regular sexual intercourse. The mean Index of Erectile Function score was 20.2 out of 25 (standard deviation [SD] 7.9), the mean Self-Esteem and Relationship score was 84.5 out of 100 (SD 9.9) and the mean Erectile Dysfunction Inventory of Treatment Satisfaction score was 82 out of 100 (SD 17.5). 13 patients (92.8%) were satisfied or very satisfied with the prosthesis.

According to the authors of the study, the results on the efficacy and safety of the prosthesis are encouraging, but further long-term and larger case studies are needed to define the clinical profile of this prosthesis.

More recently, a trial was published that included 57 patients undergoing phalloplasty surgery with a median follow-up of 16 months [2]. The trial results reported that no intraoperative complications occurred, while postoperative complications were observed in 18 of the 57 patients (32%) and included cylinder protrusion (1/57, 1.8%), infection (8/57, 14%), malposition (4/57, 7.0%), mechanical failure (4/57, 7.0%) and urinary retention (1/57, 1.8%). A total of 13 out of 57 devices (23%) were explanted, which corresponds to an explantation-free survival of 84%, 80% and 80% after 6, 12 and 18 months of follow-up, respectively.

The main limitations of this trial relate to the small sample size and the short follow-up. The rates of cylinder protrusion, malpositioning and mechanical failure are encouraging, while device infection is relatively high.

Pubmed search (<https://pubmed.ncbi.nlm.nih.gov/>); keyword: 'ZSI 475 AND phalloplasty'; search date: 18 July 2023.

Ongoing clinical trials

The search on ClinicalTrials.gov (<https://clinicaltrials.gov/>, keyword 'ZSI 475', date 18/08/2023) selected two single-arm observational studies (NCT04064671 and NCT05883553) in the process of enrolling patients.

Guidelines

Not available.

Clinical and economic evidence

Cost-effectiveness analysis

Not available.

Pubmed search (<https://pubmed.ncbi.nlm.nih.gov/>); keyword: 'ZSI 475 AND (cost[titl] OR economic[titl])'; search date: 18 July 2023.

HTA Report

Not available.

Expected benefits

The ZSI 475 FTM phalloplasty penile prostheses are the only ones that are specifically indicated for gynoandroid surgery. The technical features of these prostheses should solve the problems related to the implantation of erectile devices in the neofallus and lead to better surgical results in FtM transgender patients.

Price and cost of therapy per patient

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)	Annual requirements (No. pieces)
ZEPHYR THREE-COMPONENT PENILE PROSTHESIS (Zephyr)	8,000 + 4% VAT	8,000 + 4% VAT	4

Price and cost of therapy per patient with existing treatment alternatives

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)
No other penile prostheses suitable for phalloplasty are available	-	-

Economic and organisational impact

The purchase of the ZSI 475 FTM entails a negligible increase in expenditure given the demand for only four patients per year.

Information on reimbursability

Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM code)	DRG code	Tariff (euro)
-	-	-	-

Assessment of innovativeness (according to Regional Resolution N° 737/2022, [link](#) and N° 1244/2022 [link](#))

Innovative device (Y/N)

N

If yes, please indicate which Criteria 1, 2 and 3 are met:

Summary data

Request No.	Date of request	Applicant
286	05/04/2023	Urology - AOUC
Technology		
ZEPHYR ZSI 475 THREE-COMPONENT PENILE PROSTHESIS		

Summary data

Conclusions

Obtaining erectile function after phalloplasty in female-to-male transsexuals requires the implantation of a device. Various erectile implants exist for biological males, and numerous trials have demonstrated the success of these devices [3]. However, there are substantial differences between a native penis and a neo-phallus that explain the poorer results of these devices in female-to-male transsexuals.

Until recently, all penile prostheses were manufactured for erectile dysfunction in biological males and involved two corpora cavernosa. A few years ago, however, an inflatable 3-piece erectile device was designed specifically for female-to-male transsexuals: the Zephyr ZSI 475 FTM. This prosthesis has a realistic gland shape, a single cylinder and a specific titanium and silicone anchoring device to anchor the system to the patient's pubic bone. These specific features should solve the problems related to implanting erectile devices in the neofallus and hopefully lead to better surgical results in these patients.

The evidence supporting the ZSI 475 FTM prosthesis is still limited and the results, although encouraging, report the occurrence of complications such as infections, mechanical failure and malpositioning [1,2]. However, the request for the ZSI 475 FTM prosthesis is part of the path of assistance to people with problems related to gender identity for which the Regional Coordination Centre for Health Problems related to Gender Identity (CRIG - Regional Tuscany Region Resolution no. 907 of 6/08/2018) was established at the AOUC. The Operations Centre therefore expresses a favourable opinion.

Date of report

23/08/2023

Author of the report

Sabrina Trippoli

Referring pharmacist

Federica Romano

BIBLIOGRAPHY

1. Neuville P, Morel-Journel N, Cabelguenne D, Ruffion A, Paparel P, Terrier JE. First Outcomes of the ZSI 475 FtM, a Specific Prosthesis Designed for Phalloplasty. J Sex Med. 2019 Feb;16(2):316-322. doi: 10.1016/j.jsxm.2018.11.013. Epub 2019 Jan 9. PMID: 30638882.
2. Verla W, Goedertier W, Lumen N, Spinoit AF, Waterloos M, Waterschoot M, Buncamper M, Monstrey S, Hoebeke P. Implantation of the ZSI 475 FTM Erectile Device After Phalloplasty: A Prospective Analysis of Surgical Outcomes. J Sex Med. 2021 Mar;18(3):615-622. doi: 10.1016/j.jsxm.2020.11.008. Epub 2021 Jan 19. PMID: 33349569.
3. Preto M, Blecher G, Timpano M, Gontero P, Falcone M. The Frontier of Penile Implants in Phalloplasty: Is the ZSI 475 FTM what we have been waiting for? Int J Impot Res. 2020 Nov;33(7):779-783. doi: 10.1038/s41443-020-00396-2. Epub 2021 Jan 22. PMID: 33483605; PMCID: PMC8589655.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Drafted by the Operations Centre, Regional Decree No. 17610 of 7 September 2022.

For further information write to centro.operativo.htart@regione.toscana.it