**RAPID HTA REPORT**

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-Studio di Romano et al.:

Chronic wound is characterized by slow healing time, persistence, and abnormal healing progress. Therefore, serious complications can lead at worst to the tissue removal. In this scenario, there is an urgent need for an ideal dressing capable of high absorbency, moisture retention and antimicrobial properties. Herein we investigate the technical properties of a novel advanced non-woven triple layer gauze imbued with a cream containing Rigenase, an aqueous extract of *Triticum vulgare* used for the treatment of skin injuries. To assess the applicability of this system we analyzed the dressing properties by wettability, dehydration, absorbency, Water Vapor Transmission Rate (WVTR), lateral diffusion and microbiological tests. The dressing showed an exudate absorption up to 50%. It created a moist environment allowing a proper gaseous exchange as attested by the WVTR and a controlled dehydration rate. The results candidate the new dressing as an ideal medical device for the treatment of the chronic wound repairing process. It acts as a mechanical barrier providing a good management of the bacterial load and proper absorption of abundant wound exudate. Finally, its vertical transmission minimizes horizontal diffusion and side effects on perilesional skin as maceration and bacterial infection.

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Objectives: Compare the efficacy and tolerability of Connettivina® Bio Plus (Group A) gauze and cream, and Fitostimoline® Plus (Group B) gauze and cream for the treatment of acute superficial skin lesions. Design: Single-center, parallel, randomized trial. A block randomization method was used. Setting: University of Salerno—AOU San Giovanni di Dio e Ruggi d'Aragona.

Participants: **Sixty patients** were enrolled. All patients fulfilled the study requirements. Intervention: One application of the study drugs every 24 h, and **a six-week observation period**. Main outcome measures: Efficacy and tolerability of the study drugs.

Results: In total, 60 patients (**Group A, n = 30; Group B, n = 30**) were randomized; mean age was  $58.5 \pm 15.8$  years. All patients were included in the outcome analysis. Total wound healing was achieved in 17 patients undergoing treatment with Connettivina® Bio Plus and 28 patients undergoing treatment with Fitostimoline® Plus. The greater effectiveness of the latter was significant ( $p = 0.00104$ ). In Group B, a significantly greater degree of effectiveness was observed in reducing the fibrin in the wound bed ( $p = 0.04746$ ). Complications or unexpected events were not observed. Conclusions: Both Connettivina® Bio Plus and Fitostimoline® Plus are secure and effective for treating acute superficial skin lesions. Fitostimoline® Plus was more effective than Connettivina® Bio Plus in wound healing of acute superficial skin lesions, especially if fibrin had been observed in the wound bed.

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RAPID HTA REPORT		
Request No.	Date of request	Applicant
302	18-10-2023	UOC DERMATOLOGY, SIENA
Type of report		
New report		X
Updating a previous report		
If updating, please state the reason:		

General technology data					
<b>Trade name</b>					
PHYTOSTIMOLINE PLUS ADVANCE GAUZE					
<b>Generic name</b>					
Hyaluronated gauze impregnated with Regenerase					
<b>Manufacturer name</b>					
DAMOR Pharmaceuticals					
<b>Supplier name</b>					
Farmaceutici Damor SpA					
<b>RDM</b>	<b>REF</b>				
2115989	076701				
<b>Type</b>	<b>CE marking (date)</b>	<b>Risk class</b>	<b>FDA approval</b>		
1	6.2.2020 MED-003868-00	IIb	No		
<b>CND</b>					
M04049: DRESSINGS FOR WOUNDS, SORES AND ULCERS - OTHER					
<b>Clinical problem and rationale for the request</b>					
Dressing of various types of wounds and ulcers					
<b>Indications for use</b>					
Dressing of various types of wounds and ulcers					
<b>Target patient</b>					
Patients with a) chronic decubitus ulcers of venous or mixed vascular origin; b) surgical dehiscence; c) inflammatory ulcers; d) first degree burns; e) second degree burns; f) post-traumatic injuries; g) abrasions.					

### General technology data

#### Main competitors

It is emphasised that the applicant proposes to use DM in 2 consecutive phases: in the first phase, advanced dressings with antiseptic properties are used, followed by hyaluronate gauze

### Technological details

#### Description

Fitostimoline Plus gauze is a medical device for dermatological use consisting of single-dose, PET gauze pads, 10x10 cm in size, impregnated with water-dispersible cream containing Rigenase® and Polyhexanide, individually packed in sealed paper/aluminium/PE bags. The product, in its final packaging, undergoes a terminal  $\gamma$ -ray sterilisation process. The medical device Fitostimoline® Plus Advance gauze applied to the injured skin forms a protective barrier against the external environment contributing to the control of the wound microenvironment and thus creating favourable conditions for a rapid and correct re-epithelialising action on the skin. The product is therefore indicated for the treatment of wounds of different aetiologies. The special triple PET layer support allows the management of moderate/intensive exudates. Rigenase® is a specific wheat extract. The presence of polyhexanide reduces the risk of bacterial contamination.

#### Elements of Innovation

None in particular

### Clinical and economic evidence

#### Clinical studies

There are three main studies supporting this DM [1-3].

The study by Ricci et al [1], on the one hand, was published in an Italian journal not listed in PubMed and, on the other hand, only included 20 patients in a single-arm trial, i.e. conducted in the absence of a control group. The analysis towards a control group was surrogated by an 'after vs. before' comparison for each patient compared with himself.

The study by Romano et al. [2] is pre-clinical in vitro and therefore did not enrol any patients. It reports the main parameters measured in vitro in the pharmacological development phase of DM.

The study by Russo et al. [3] is the most valid one, but it still only enrolled 30 patients treated with Phytostimulines for 6 weeks.

#### Ongoing clinical trials

The clinicaltrials.gov database as of 6 November 2023 reported only one ongoing study in patients with anal fissure; no results available.

#### Guidelines

-no

#### Cost-effectiveness analysis

-no

#### HTA Report

-none

#### Expected benefits

**Clinical and economic evidence**

-unquantified/unquantifiable

**Price and cost of therapy per patient**

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)	Annual requirements (No. pieces)
PHYTOSTIMOLINE PLUS ADVANCE GAUZE (DAMOR)	29.60 € in a pack of 5 gauze (corresponding to approximately 6 € per gauze)	It depends on the duration of therapy with this DM; the study by Russo et al (3) who used a duration of 6 weeks.DM	about 450 gauze

**Price and cost of therapy per patient with existing treatment alternatives**

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)
The applicant proposes an add-on use, because an active therapeutic alternative cannot be identified.		

**Economic and organisational impact**

Approximately €3,000 per year limited to the applicant structure

**Reimbursability information**

Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM Code )	DRG code	Tariff (euro)
A variety of clinical conditions are candidates for therapy with this DM		Unidentifiable	

**Assessment of innovativeness (according to Regional Resolution N° 737/2022, [link](#) and N° 1244/2022 [link](#) )**

Innovative device (Y/N)	N
If yes, please indicate which Criteria 1, 2 and 3 are met:	

**Summary data**

Request number	Date of request	Applicant
302	18-10-2023	AOUS Dermatology
Technology		
Device for use as a dressing		

### Summary data

The application for Phytostimoline PLUS ADVANCE is essentially an off-label use of DM since the dosage proposed in the two available clinical trials provides therapy for 6 weeks (thus different from the variable duration dosage that the applicant intends to use). In principle, this does not constitute an absolute reason against a positive recommendation. However, it does pose the problem that any long-term therapy (which does not seem to correspond to the applicant's intentions in any case) would inevitably include a period of administration after the patient has been discharged, in which the financial burden of the purchase falls on the patient. A possible approval could be justified by the modest economic impact this DM would have. However, the framework for the use of DM appears ill-defined and also different from the clinical studies supporting DM itself.

In conclusion, given that the evidence is scanty and in any case refers to the use of Phytostimoline for 6 weeks, the CO expresses an unfavourable opinion considering the possible evidence in favour of DM to be insufficient. In this decision, the unit price of Phytostimoline PLUS, which appears considerably higher than that of other medicated gauze, plays a major role.

### Date of report

14 November 2023

Andrea Messori

### Referring pharmacist

Donata Iozzi

### BIBLIOGRAPHY

1. Ricci E, Pittarello M. Clinical evaluation of the efficacy and tolerability of a new Rigenase® and polyhexanide-based dressing in the treatment of chronic skin lesions **Italian Journal of Wound Care**, 2022, № 3, <https://doi.org/10.4081/ijwc.2022.92>
2. Romano E, Campagnuolo C, Palladino R, Schiavo G, Maglione B, Luceri C, Mennini N. Technical Evaluation of a New Medical Device Based on Rigenase in the Treatment of Chronic Skin Lesions. **Bioengineering** (Basel). 2023 Aug 29;10(9):1022. doi: 10.3390/bioengineering10091022. PMID: 37760124; PMCID: PMC10526047.
3. Russo R, Carrizzo A, Barbato A, Rasile BR, Pentangelo P, Ceccaroni A, Marra C, Alfano C, Losco L. Clinical Evaluation of the Efficacy and Tolerability of Rigenase® and Polyhexanide (Fitostimoline® Plus) vs. Hyaluronic Acid and Silver Sulfadiazine (Connettivina® Bio Plus) for the Treatment of Acute Skin Wounds: A Randomized Trial. **J Clin Med**. 2022 Apr 29;11(9):2518. doi: 10.3390/jcm11092518. PMID: 35566643; PMCID: PMC9105357.

## APPENDIX: ABSTRACTS OF THE ABOVE REFERENCES

-Study by Ricci et al.

This article discusses the case control of a dressing composed of Rigenase® (plant-derived stimulin) with added polyhexanide, on a new triple-layer PET carrier in the treatment of chronic skin lesions of the lower limb. The aim of the study is the evaluation in terms of repair and control of bacterial bioburden. **Twenty patients were** enrolled with a run-in period of 6 weeks with gold standard therapy, followed by a similar period of treatment with the investigational dressing. The end points were the evolution of the area, measured with an artificial intelligence-based system as primary. Secondary end points: level of infection/colonisation, WBP and pain, all collected using validated scales. The data collected with regard to area showed a statistically significant difference ( $p < 0.0001$ ). Significant differences were also evident in terms of signs of infection ( $p < 0.001$ ), WBP score ( $p < 0.001$ ) and pain with NRS ( $p < 0.0001$ ). No adverse events or side effects were observed. The authors conclude that the dressing demonstrates a high capacity to stimulate the resumption of reparative phenomena in lesions with a poor tendency to spontaneous repair.

-Study by Romano et al:

Chronic wound is characterised by slow healing time, persistence, and abnormal healing progress. Therefore, serious complications can lead at worst to the tissue removal. In this scenario, there is an urgent need for an ideal dressing capable of high absorbency, moisture retention and antimicrobial properties. Herein we investigate the technical properties of a novel advanced non-woven triple layer gauze imbibed with a cream containing Rigenase, an aqueous extract of *Triticum vulgare* used for the treatment of skin injuries. To assess the applicability of this system we analysed the dressing properties by wettability, dehydration, absorbency, Water Vapor Transmission Rate (WVTR), lateral diffusion and microbiological tests. The dressing showed an exudate absorption up to 50%. It created a moist environment allowing a proper gaseous exchange as attested by the WVTR and a controlled dehydration rate. The results candidate the new dressing as an ideal medical device for the treatment of the chronic wound repairing process. It acts as a mechanical barrier providing a good management of the bacterial load and proper absorption of abundant wound exudate. Finally, its vertical transmission minimises horizontal diffusion and side effects on perilesional skin as maceration and bacterial infection.

-Study by Russo et al.

Objectives: Compare the efficacy and tolerability of Connettivina® Bio Plus (Group A) gauze and cream, and Fitostimoline® Plus (Group B) gauze and cream for the treatment of acute superficial skin lesions. Design: Single-center, parallel, randomised trial. A block randomisation method was used. Setting: University of Salerno-AOU San Giovanni di Dio e Ruggi d'Aragona.

Participants: **Sixty patients** were enrolled. All patients fulfilled the study requirements. Intervention: One application of the study drugs every 24 h, and **a six-week observation period**. Main outcome measures: Efficacy and tolerability of the study drugs.

Results: In total, 60 patients (**Group A, n = 30; Group B, n = 30**) were randomised; mean age was  $58.5 \pm 15.8$  years. All patients were included in the outcome analysis. Total wound healing was achieved in 17 patients undergoing treatment with Connettivina® Bio Plus and 28 patients undergoing treatment with Fitostimoline® Plus. The greater effectiveness of the latter was significant ( $p = 0.00104$ ). In Group B, a significantly greater degree of effectiveness was observed in reducing the fibrin in the wound bed ( $p = 0.04746$ ). Complications or unexpected events were not observed. Conclusions: Both Connettivina® Bio Plus and Fitostimoline® Plus are secure and effective for treating acute superficial skin lesions. Fitostimoline® Plus was more effective than Connettivina® Bio Plus in wound healing of acute superficial skin lesions, especially if fibrin had been observed in the wound bed.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Reporting by the Operations Centre, Regional Decree No. 17610 of 7 September 2022.

For further information write to [centro.operativo.htart@regione.toscana.it](mailto:centro.operativo.htart@regione.toscana.it)