**RAPID HTA REPORT**

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
| RAPID HTA REPORT                      |                 |                                      |
|---------------------------------------|-----------------|--------------------------------------|
| Request No.                           | Date of request | Applicant                            |
| 296                                   | 11/07/2023      | UOC Interventional Cardiology - FTGM |
| Type of report                        |                 |                                      |
| New report                            |                 | Yes                                  |
| Updating a previous report            |                 | No                                   |
| If updating, please state the reason: |                 |                                      |

| General technology data   |                                   |   |                     |
|---|-----------------------------------|---|---------------------|
| <b>Trade name</b>   |                                   |  |                     |
| TricValve Bicaval Valve System SVC25, SVC29, IVC31, IVC35   |                                   |   |                     |
| <b>Generic name</b>   |                                   |   |                     |
| Bicaval transcatheter prosthesis for the treatment of tricuspid insufficiency   |                                   |   |                     |
| <b>Manufacturer name</b>  |                                   |   |                     |
| P+F PRODUCTS + FEATURES GMBH  |                                   |   |                     |
| <b>Supplier name</b>  |                                   |   |                     |
| Innova HTS srl  |                                   |   |                     |
| <b>RDM</b>  | <b>REF</b>                        |   |                     |
| 2113482   | 1101020025                        |   |                     |
| 2113543   | 1101020029                        |   |                     |
| 2113544   | 1101020031                        |   |                     |
| 2113545   | 1101020035                        |   |                     |
| <b>Type</b>   | <b>CE marking (date)</b>          | <b>Risk class</b>   | <b>FDA approval</b> |
| 1   | 1434-MDD-213/2021<br>(14/05/2021) | III   | No                  |
| <b>CND</b>  |                                   |   |                     |
| P0799 (VASCULAR AND CARDIAC PROSTHESES-OTHER)   |                                   |   |                     |
| <b>Clinical problem and rationale for the request</b>   |                                   |   |                     |
| <p>Severe symptomatic tricuspid valve regurgitation (RT) is associated with inadequate clinical outcomes in the short to medium term, representing the major cause of moderate to severe valve disease in developed countries. Medical treatment and valve surgery are the currently accepted therapies to treat severe symptomatic tricuspid valve (VT) regurgitation. However, although volumes of tricuspid valve surgery are steadily increasing, the hospital mortality rate remaining at 9% suggests that the surgical solution is offered too late in the disease course [2-3].</p> <p>Several percutaneous devices have emerged in the last decade for the treatment of RT, with promising results in feasibility studies and ongoing clinical trials. However, given the anatomical complexity and peri-procedural imaging challenges, many patients are still inoperable with percutaneous treatments (edge-to-edge repair, annuloplasty, orthotopic valve replacement) [4-5]. Caval valve implantation (CAVI) initially emerged as a therapeutic alternative for patients referred to as "with no other options" for the treatment of symptomatic and severe RT; these patients usually also suffer from liver congestion and heart failure [6].</p> |                                   |   |                     |

### General technology data

Given the emerging challenges and applicability of many emerging percutaneous therapies to the broader cohort of candidates with severe valvular regurgitation, CAVI-implanted procedures have some inherent advantages due to the relative simplicity of the caval anatomy compared to VT and right heart [7] and the anatomical and peri-procedural imaging requirements for CAVI are less stringent.

Therefore, when possible, it is preferred not to intervene directly on the VT but more simply on the caval area.

### Indications for use

The TricValve device falls into the category of heterotopic caval valve implants (CAVI). The bioprosthesis was developed for the treatment of the caval reflux that accompanies cases of severe tricuspid regurgitation in high-risk open-heart surgery. The bioprosthesis was designed to treat severe tricuspid regurgitation without removing the defective tricuspid valve.

### Target patient

Patients with severe symptomatic tricuspid valve regurgitation at high surgical risk or inoperable

### Main competitors

The TricValve is currently the only CE-marked CAVI device on the market that allows treatment of the target population without the need for tricuspid valve repair. If repair techniques cannot be used, the pathological tricuspid valve must be treated with surgical replacement techniques (associated with poor outcomes) or transcatheter replacement (not yet on the market).

### Technological details

#### Description

The bioprosthesis is available in two different diameters for each model (VCS and VCI) and is specifically designed to fit the anatomical features of the superior and inferior vena cava. The bioprosthesis consists of a metal tubular nitinol structure, which is self-expanding and radiopaque, and is equipped with three sutured bovine pericardium valve flaps complemented by a polyester collar to prevent para-valvular leakage. The flaps of the bioprosthesis are treated with an anti-calcification product and subjected to chemical dehydration.

#### Elements of Innovation

CAVI procedures have some inherent advantages due to the relative simplicity of the caval anatomy compared to the tricuspid valve and right heart. The anatomical (and peri-procedural imaging) requirements for CAVI are less stringent. Therefore, when possible, it is preferred not to operate directly on the tricuspid valve but more simply on the caval area.

### Clinical and economic evidence

#### Clinical studies

The PubMed search with the keyword 'TricValve' returned 17 results from which 2 studies were selected:

#### **Interventional Treatment of Severe Tricuspid Regurgitation Early Clinical Experience in a Multicenter, Observational, First-in-Man Study [8]**

25 patients (females: 52%; age,  $73.9 \pm 7.6$  years; EuroScore II,  $18.2 \pm 12.9$ ; STS  $14.0 \pm 12.7$ ) with massive TR and symptoms of right heart failure were selected for this new therapeutic approach within the framework of a compassionate clinical use programme. At baseline, 72% (n = 18) and 28% (n = 7) of the patients were in NYHA class IV and III, respectively. Functional TR was present in 96% (n = 24) and structural TR (rheumatic) in 4% (n = 1) of patients. The total number of valves implanted within this study was 31. Patients were treated with implantation of a single valve (IVC only; n=19; 76.0%) or double BiCAVI (IVC and SVC; n=6; 24.0%). In these patients, balloon

**Clinical and economic evidence**

expandable valves (Sapien XT or Sapien 3: n=17; 78.3%) or self-expanding valves (TricValve: n=7; 21.7%; DirectFlow: n=1; 4%) were used. The 30-day mortality was 12% (n=3). Death was caused by progressive multi-organ failure and complications due to sepsis after 6, 14 and 24 days after surgery. In-hospital mortality was 24% (n=6). Two patients developed pneumonia and died from respiratory complications 34 and 41 days after implantation and one patient died from multi-organ failure 45 days after implantation. The mean follow-up was 316±453 days (14-1540 days). In all patients discharged alive (n=19), transthoracic echocardiography during follow-up confirmed overall function of all implanted devices. The longest follow-up (51 days) was one Bi CAVI case with documented intact device function. At discharge, the patients' symptoms improved in 84.2% (n=16) and there was also improvement in one or more NYHA classes. At discharge 9 patients had NYHA class III or IV (47.3%) and 52.7% improved by one or two classes.

**6-Month Outcomes of the TricValve System in Patients With Tricuspid Regurgitation: The TRICUS EURO Study [9].**

35 patients (mean age 76 ± 6.8 years, 83% women) were treated with the TricValve system. At baseline, all patients had NYHA functional class III or IV. At 30 days, procedural success was 94%, with no procedural death or conversion to surgery. At 6-month follow-up, a significant increase in QOL (Kansas City Cardiomyopathy Questionnaire scores at baseline and at 6 months, 42.01 ± 22.3 and 59.7 ± 23.6, respectively; p = 0.004) was observed, correlating with a significant improvement in NYHA class, with 79.4% of patients in NYHA class I or II at 6 months (p = 0.0006). At 6 months, the all-cause mortality and heart failure hospitalisation rates were 8.5% and 20%, respectively.

**Ongoing clinical trials**

The search on clinicaltrials.gov with the keyword 'TricValve' revealed 4 studies:

- NCT05820516: international registry of 100 participants with severe RT undergoing TricValve implantation from Italian and European centres (not yet recruited).
- NCT04141137: Single-arm interventional study of 35 patients with severe RT and significant reflux in the superior and/or inferior vena cava who underwent TricValve implantation. The study was completed and the results were published at 6 months [8].
- NCT03723239: single-arm interventional study of 9 patients with severe RT and significant reflux in the superior and/or inferior vena cava who underwent TricValve implantation. The study was completed but no results appear to have been published.
- NCT05114850: register of 450 participants treated with TricValve implant (under recruitment).

**Guidelines**

-

**Cost-effectiveness analysis**

-

**HTA Report**

-Report Ontario:

<https://database.inahta.org/article/20974>

Based on guidance from the Ontario Health Technology Advisory Committee, the HTA report recommends public funding of transcatheter valve implantation in adult patients with degenerated mitral bioprosthetic valves (TMViV) or tricuspid valves (TTViV), which are considered inoperable or at high risk for surgery. The report concludes by estimating that public funding of TMViV and TTViV in Ontario would result in savings of \$0.33 million in 5 years.

**Expected benefits**

Reduction of reflux in the superior and/or inferior vena cava related to improved RT.

**Price and cost of therapy per patient**

| Product (Manufacturer)  | Unit price (euro) | Therapy cost per patient (euro) | Annual requirements (No. pieces) |
|---|-------------------|---------------------------------|----------------------------------|
| TricValve Bicaval Valve System (P+F PRODUCTS + FEATURES GMBH) | 24.000,00         | 24.000,00                       | 5                                |

**Price and cost of therapy per patient with existing treatment alternatives**

| Product (Manufacturer) | Unit price (euro) | Therapy cost per patient (euro) |
|------------------------|-------------------|---------------------------------|
| Surgery                |                   |                                 |

**Economic and organisational impact**

|  |
|--|
|  |
|--|

**Reimbursability information**

| Main diagnosis (ICD9-CM code)                                | Intervention (ICD9-CM Code )                           | DRG code | Tariff (euro) |
|--|--|----------|---------------|
| 397.0 (Tricuspid valve disease)                              | 35.27 (Tricuspid valve replacement with bioprosthesis) | 104      | 24.115,00     |
| 424.2 (Tricuspid valve disorders specified as non-rheumatic) | 38.91 (Arterial catheterisation)                       |          |               |

**Assessment of innovativeness (according to Regional Resolution N° 737/2022, [link](#) and N° 1244/2022 [link](#) )**

|  |   |
|--|---|
| Innovative device (Y/N)                                    | N |
| If yes, please indicate which Criteria 1, 2 and 3 are met: |   |

**Summary data**

| Request number  | Date of request | Applicant                            |
|---|-----------------|--------------------------------------|
| 296   | 11/07/2023      | UOC Interventional Cardiology - FTGM |
| Technology  |                 |                                      |
| TricValve Bicaval Valve System SVC25, SVC29, IVC31, IVC35 (Bicaval transcatheter prosthesis for the treatment of tricuspid insufficiency) |                 |                                      |

### Summary data

Due to the high operative mortality of isolated tricuspid valve surgery in the case of late operations [2-3], it is probably fair to predict that the percutaneous approach will be the treatment of choice in high-risk and/or inoperable patients requiring surgery on this valve.

Considering the great therapeutic need for an interventional option for inoperable patients with severe tricuspid regurgitation, percutaneous valve implantation might be a feasible alternative, in addition to the now more widespread percutaneous repair.

Compared to other devices with a tricuspid indication, procedures performed with caval valve implantation (CAVI) have some inherent advantages due to the relative simplicity of the caval anatomy compared to the tricuspid valve and right heart [3] and the anatomical (and peri-procedural imaging) requirements for CAVI are less stringent, the TricValve device is unique in having CE marking with this indication.

The TricValve system is suitable for a target population at high surgical risk or inoperable with tricuspid valve insufficiency and caval reflux.

### Date of report

02/10/23

### Author(s) of the report

Andrea Messori, Maria Rita Romeo

### Referring pharmacist

Stefania Alduini

## BIBLIOGRAPHY

1. Romeo MR, Nardone A. Rapid Report HTA System TricValve (preprint). Open Science Framework, DOI 10.17605/OSF.IO/NSV9W, url <https://osf.io/u2mt9>, accessed 21 August 2023.
2. Topilsky Y, Maltais S, Medina Inojosa J, Oguz D, Michelena H, Maalouf J, Mahoney DW, Enriquez-Sarano M. Burden of Tricuspid Regurgitation in Patients Diagnosed in the Community Setting. *JACC Cardiovasc Imaging*. 2019 Mar;12(3):433-442. doi: 10.1016/j.jcmg.2018.06.014. Epub 2018 Aug 15. PMID: 30121261.
3. Stewart WJ. Tricuspid Regurgitation: Moderate (or More) Is a Burden for the Right Ventricle. *JACC Cardiovasc Imaging*. 2019 Mar;12(3):443-445. doi: 10.1016/j.jcmg.2018.07.018. Epub 2019 Jan 16. PMID: 30660552.
4. Zack CJ, Fender EA, Chandrashekar P, Reddy YNV, Bennett CE, Stulak JM, Miller VM, Nishimura RA. National Trends and Outcomes in Isolated Tricuspid Valve Surgery. *J Am Coll Cardiol*. 2017 Dec 19;70(24):2953-2960. doi: 10.1016/j.jacc.2017.10.039. PMID: 29241483.
5. Kim JB, Jung SH, Choo SJ, Chung CH, Lee JW. Surgical outcomes of severe tricuspid regurgitation: predictors of adverse clinical outcomes. *Heart*. 2013 Feb;99(3):181-7. doi: 10.1136/heartjnl-2012-302856. Epub 2012 Oct 4. PMID: 23038792.

6. Asmarats L, Puri R, Latib A, Navia JL, Rodés-Cabau J. Transcatheter Tricuspid Valve Interventions: Landscape, Challenges, and Future Directions. *J Am Coll Cardiol*. 2018 Jun 26;71(25):2935-2956. doi: 10.1016/j.jacc.2018.04.031. PMID: 29929618.
7. Lauten A, Figulla HR, Willich C, Jung C, Rademacher W, Schubert H, Ferrari M. Heterotopic valve replacement as an interventional approach to tricuspid regurgitation. *J Am Coll Cardiol*. 2010 Feb 2;55(5):499-500. doi: 10.1016/j.jacc.2009.09.034. PMID: 20117467.
8. Lauten A, Figulla HR, Unbehaun A, Fam N, Schofer J, Doenst T, Hausleiter J, Franz M, Jung C, Dreger H, Leistner D, Alushi B, Stundl A, Landmesser U, Falk V, Stangl K, Laule M. Interventional Treatment of Severe Tricuspid Regurgitation: Early Clinical Experience in a Multicenter, Observational, First-in-Man Study. *Circ Cardiovasc Interv*. 2018 Feb;11(2):e006061. doi: 10.1161/CIRCINTERVENTIONS.117.006061. PMID: 29445001.
9. Estévez-Loureiro R, Sánchez-Recalde A, Amat-Santos IJ, Cruz-González I, Baz JA, Pascual I, Mascherbauer J, Abdul-Jawad Altisent O, Nombela-Franco L, Pan M, Trillo R, Moreno R, Delle Karth G, Salido-Tahoces L, Santos-Martinez S, Núñez JC, Moris C, Goliasch G, Jimenez-Quevedo P, Ojeda S, Cid-Álvarez B, Santiago-Vacas E, Jimenez-Valero S, Serrador A, Martín-Moreiras J, Strouhal A, Hengstenberg C, Zamorano JL, Puri R, Íñiguez-Romo A. 6-Month Outcomes of the TricValve System in Patients With Tricuspid Regurgitation: The TRICUS EURO Study. *JACC Cardiovasc Interv*. 2022 Jul 11;15(13):1366-1377. doi: 10.1016/j.jcin.2022.05.022. Epub 2022 May 17. PMID: 35583363.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Reporting by the Operations Centre, Regional Decree No. 17610 of 7 September 2022.

For further information write to [centro.operativo.htart@regione.toscana.it](mailto:centro.operativo.htart@regione.toscana.it)