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4. Oh D, Lee JH, Song TJ, Song KB, Hwang DW, Kim JH, Park DH, Lee SS, Seo DW, Lee SK, Kim MH. Clinical outcomes of EUS-guided transluminal drainage with a novel lumen-apposing metal stent for postoperative pancreatic fluid collection after pancreatic surgery. *Gastrointest Endosc.* 2022 Apr;95(4):735-746. doi: 10.1016/j.gie.2021.12.015. Epub 2021 Dec 28. PMID: 34971669.
5. Mangiavillano B, Auriemma F, Lamonaca L, Repici A. A Novel Lumen-Apposing Metal Stent with an Electrocautery Tip for Different Indications: Initial Experience in a Referral Center. *Dig Dis.* 2022;40(4):526-529. doi: 10.1159/000518539. Epub 2021 Aug 4. PMID: 34348270; PMCID: PMC8450863.
6. Teoh AYB, Bapaye A, Lakhtakia S, Ratanachu T, Reknimitr R, Chan SM, Choi HJ, Gadhikar HP, Kongkam P, Korrapati SK, Lee YN, Medarapalem J, Ridditid W, Moon JH. Prospective multicenter international study on the outcomes of a newly developed self-approximating lumen-apposing metallic stent for drainage of pancreatic fluid collections and endoscopic necrosectomy. *Dig Endosc.* 2020 Mar;32(3):391-398. doi: 10.1111/den.13494. Epub 2019 Aug 30. PMID: 31343773.
7. Song TJ, Lee SS, Moon JH, Choi HJ, Cho CM, Lee KH, Park SW, Kim SH, Lee SO, Lee YN, Lee JK. Efficacy of a novel lumen-apposing metal stent for the treatment of symptomatic pancreatic pseudocysts (with video). *Gastrointest Endosc.* 2019 Sep;90(3):507-513. doi: 10.1016/j.gie.2019.05.033. Epub 2019 May 30. PMID: 31152738.

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**RAPID HTA REPORT**

Request No.	Date of request	Applicant
303	30/05/2023	UOC Digestive Endoscopy - Livorno Hospital
Type of report		
New report		X
Updating a previous report		
If updating, please state the reason:		

General technology data

Trade name				
Niti-S Hot SPAXUS™ Stents				
Generic name				
Stents with electrocautery delivery system				
Manufacturer name				
Taewoong Medical				
Supplier name				
Euromedical				
RDM				REF
1875371				HSS1602FW
Type	CE marking (date)	Risk class	FDA approval	
1	KR19/81826260 (22/10/2019)	IIb	Yes	

CND

P0502 (Biliary and pancreatic prostheses)

Clinical problem and rationale for the request

The rationale for the request is based on an economic aspect. A stent that can be superimposed on the Niti-S Hott Spaxus, the Hot Axios stent (previously evaluated, link: https://www301.regione.toscana.it/bancadati/atti/Contenuto.xml?id=5262373&nomeFile=Decreto_n.12735_del_13-08-2020-Allegato-6), is currently in use. The Hot Axios stent has a unit price of €3,000, while the Hot Spaxus stent has a unit price of €2,500. The approval of the Niti-S Hot Spaxus would result in a saving per procedure of €1,000 considering that two stents are used per patient. The Niti-S Hot Spaxus stent is requested as a supplement to the Hot Axios and not as a replacement for it.

Indications for use

The Niti-S Hot SPAXUS™ stent with electrocautery delivery system [HDS-10] is designed for the drainage of pancreatic pseudocysts or gallbladder pseudocysts with a transgastric or transduodenal approach.

Target patient

Patients suffering from:

- complications of acute or chronic pancreatitis;
- stenosis of the biliary tract with obstructions that cannot be overcome by retrograde drainage;
- cholecystitis that cannot be treated surgically.

Main competitors

Other lumen-apposing metal stents (LAMS) such as the Hot Axios stent, plastic stents and covered metal stents.

Technological details

Description

The stent is made of Nitinol wire. It is a flexible, fine-mesh tubular stent with radiopaque markers on each end and in the centre.

The electrocautery stent delivery system is designed to be connected to an electrosurgical unit and penetrate adjacent human organ tissue with electrocautery. The stent is loaded into the electrocautery stent delivery system [HDS-10] and upon deployment the body of the stent imparts a radial outward force, while both ends of the device impart a horizontal inward force, allowing transmural drainage.

The Niti-S Hot SPAXUS™ stent is completely covered with silicone and can be removed.

The Electrocautery Delivery System [HDS-10] is a delivery system with an advanced electrocautery function compatible with therapeutic echo-endoscopes and a working channel diameter of 3.7 mm or more.

The endoscopic electrocautery delivery system [HDS-10] has a usable length of 180 cm. The endoscopic type is recommended when an endoscopic approach is intended.

Elements of Innovation

The Niti-S Hot Spaxus stent does not present any technological innovations compared to the Hot Axios stent, however, the Hot Axios stent has two different placement systems in the target structure, Seldinger technique and Freestyle technique, whereas the Hot Spaxus stent has only one placement system with Seldinger technique.

Clinical and economic evidence

Clinical studies

The Niti-S Hot Spaxus stent has been evaluated in several clinical studies with a single-arm observational design [1-7]. The studies have generally enrolled a small cohort of patients with short follow-up. The most relevant trial was a propensity-matching analysis [1] that compared the Hot Spaxus stent (n=132) to the Hot Axios stent (n=132) for endoscopic drainage of pancreatic fluid. Technical and clinical success rates were comparable between the two groups. Significantly more haemorrhages requiring transfusion and/or intervention occurred in the Hot Axios group than in the Hot Spaxus group (6.8% vs. 1.5%; $P = 0.03$); stent type was a significant predictor of haemorrhage in both univariate and multivariate regression analyses ($P = 0.03$ and 0.04 , respectively). Bleeding that required arterial embolisation, however, did not differ significantly between the two groups (3.0% vs. 0%; $P = 0.12$). Furthermore, Hot Axios was associated with a significantly higher rate of overall adverse events compared to the Spaxus stent (9.8% vs. 3.0%; $P = 0.04$).

Pubmed search (<https://pubmed.ncbi.nlm.nih.gov/>): keyword: axios; search date: 23/02/2024.

Ongoing clinical trials

-

Guidelines

Not available

Cost-effectiveness analysis

Not available

HTA Report

Not available

Expected benefits

Reducing expenditure

Price and cost of therapy per patient

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)	Annual requirements (No. pieces)
---------------------------	----------------------	------------------------------------	-------------------------------------

Price and cost of therapy per patient

Niti-S Hot SPAXUS™ Stents (Taewoong Medical)	2,500 excluding VAT	about 5,000	50*
---	---------------------	-------------	-----

*There are 25 patients for whom Hot Spaxus has been requested. Two stents are used for each patient.

Price and cost of therapy per patient with existing treatment alternatives

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)
HOT Axios stents (Boston Scientific)	3,000 excluding VAT	6,000 (considering the use of two stents per patient)

Economic and organisational impact

The total annual number of pancreatic pseudocyst or gallbladder drainage procedures with a transgastric or transduodenal approach by the Livorno Hospital Endoscopy is 50. Of these 50, half would be performed with Hot Spaxus and the other half with Hot Axios. Since the price of the Niti-S Hot Spaxus stent is 500 euros less than that of the Hot Axios stent, and two stents are used for each procedure, the savings per patient from the use of Hot Spaxus is approximately 1,000 euros.

Reimursability information

Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM Code)	DRG code	Tariff (euro)
Other cholecyst disorders (575) Postcholecystectomy syndrome (576) Malignant tumours of the liver and intrahepatic bile ducts (155)	-	193 (Biliary tract operations except isolated cholecystectomy with or without exploration of the common bile duct with CC)	13.509
		194 (Biliary tract operations except isolated cholecystectomy with or without exploration of the common bile duct without CC)	8.789

Assessment of innovativeness (according to Regional Resolution N° 737/2022, [link](#) and N° 1244/2022 [link](#))

Innovative device (Y/N)	N
If yes, please indicate which Criteria 1, 2 and 3 are met:	

Summary data		
Request number	Date of request	Applicant
303	30/05/2023	UOC Digestive Endoscopy - Livorno Hospital
Technology		
Niti-S Hot SPAXUS™ stent (stent with electrocautery delivery system)		
Conclusions		
<p>The Niti-S Hot SPAXUS stent with electrocautery delivery system is designed for the drainage of pancreatic pseudocysts or gallbladder pseudocysts with a transgastric or transduodenal approach.</p> <p>The Niti-S Hot SPAXUS stent was requested as a supplement to the Hot Axios stent in use since 2018 following approval by the then regional medical device evaluation group, the GRDM.</p> <p>A propensity-matching analysis [1] reported that Hot Spaxus has a better safety profile than Hot Axios; however, further studies are needed to better define this aspect.</p> <p>In comparison with the Hot Axios stent, Hot Spaxus results in a reduction in expenditure per procedure of €1,000.</p> <p>Hot Axios and Hot Spaxus have different introduction systems into the target structure, and given the complexity for pancreatic and biliary pseudocyst drainage procedures with a transgastric or duodenal approach, it is considered necessary to be able to use both types of stent.</p> <p>Based on these considerations, the CO expresses a favourable opinion.</p>		
Date of report		
22/02/2024		
Author(s) of the report		
Sabrina Trippoli		
Referring pharmacist		
Valentina Berni		

BIBLIOGRAPHY

1. Mangiavillano B, Lakhtakia S, Samanta J, Auriemma F, Vargas-Madrigal J, Arcidiacono PG, Barbera C, Ashhab H, Song TJ, Pham KD, Teoh AYB, Moon JH, Crinò SF, Kongkam P, Aragona G, De Lusong MA, Dhar J, Ofosu A, Ventra A, Paduano D, Franchellucci G, Repici A, Larghi A, Facciorusso A; PFC LAMS study group. Lumen-apposing metal stents for the treatment of pancreatic and peripancreatic fluid collections and bleeding risk: a propensity matched study. *Endoscopy*. 2024 Jan 18. doi: 10.1055/a-2219-3179. Epub ahead of print. PMID: 38237633.
2. Mangiavillano B, Moon JH, Facciorusso A, Vargas-Madrigal J, Di Matteo F, Rizzatti G, De Luca L, Forti E, Mutignani M, Al-Lehibi A, Paduano D, Bulajic M, Decembrino F, Auriemma F, Franchellucci G, De Marco A, Gentile C, Shin IS, Rea R, Massidda M, Calabrese F, Mirante VG, Ofosu A, Crinò SF, Hassan C, Repici A, Larghi A. Endoscopic ultrasound-guided gallbladder drainage as a first approach for jaundice palliation in unresectable malignant distal biliary obstruction: Prospective study. *Dig Endosc*. 2023 May 30. doi: 10.1111/den.14606. Epub ahead of print. PMID: 37253185.
3. Manta R, Torrisi S, Castellani D, Germani U, Zito FP, Morelli O, Brunori PM, Zullo A, Francesco V. Endoscopic ultrasonography-guided biliary drainage for malignant jaundice by using electrocautery-enhanced lumen-apposing metal stents (Hot-SPAXUS stent). *Ann Gastroenterol*. 2022 Nov-Dec;35(6):663-667. doi: 10.20524/aog.2022.0743. Epub 2022 Oct 3. PMID: 36406962; PMCID: PMC9648518.
4. Oh D, Lee JH, Song TJ, Song KB, Hwang DW, Kim JH, Park DH, Lee SS, Seo DW, Lee SK, Kim MH. Clinical outcomes of EUS-guided transluminal drainage with a novel lumen-apposing metal stent for postoperative pancreatic

- fluid collection after pancreatic surgery. *Gastrointest Endosc.* 2022 Apr;95(4):735-746. doi: 10.1016/j.gie.2021.12.015. Epub 2021 Dec 28. PMID: 34971669.
5. Mangiavillano B, Auriemma F, Lamonaca L, Repici A. A Novel Lumen-Apposing Metal Stent with an Electrocautery Tip for Different Indications: Initial Experience in a Referral Center. *Dig Dis.* 2022;40(4):526-529. doi: 10.1159/000518539. Epub 2021 Aug 4. PMID: 34348270; PMCID: PMC8450863.
 6. Teoh AYB, Bapaye A, Lakhtakia S, Ratanachu T, Reknimitr R, Chan SM, Choi HJ, Gadhikar HP, Kongkam P, Korrapati SK, Lee YN, Medarapalem J, Ridditid W, Moon JH. Prospective multicenter international study on the outcomes of a newly developed self-approximating lumen-apposing metallic stent for drainage of pancreatic fluid collections and endoscopic necrosectomy. *Dig Endosc.* 2020 Mar;32(3):391-398. doi: 10.1111/den.13494. Epub 2019 Aug 30. PMID: 31343773.
 7. Song TJ, Lee SS, Moon JH, Choi HJ, Cho CM, Lee KH, Park SW, Kim SH, Lee SO, Lee YN, Lee JK. Efficacy of a novel lumen-apposing metal stent for the treatment of symptomatic pancreatic pseudocysts (with video). *Gastrointest Endosc.* 2019 Sep;90(3):507-513. doi: 10.1016/j.gie.2019.05.033. Epub 2019 May 30. PMID: 31152738.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Reporting by the Operations Centre, Regional Decree No. 17610 of 7 September 2022.

For further information write to centro.operativo.htart@regione.toscana.it