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RAPID HTA REPORT		
Request No.	Date of request	Applicant
324	30/05/2024	Cardiac Surgery - AOUC
Type of report		
New report		X
Updating a previous report		
If updating, please state the reason:		

General technology data			
<b>Trade name</b>			
IMPELLA 5.5 CATHETER KIT WITH SMART ASSIST			
<b>Generic name</b>			
Minimally invasive left ventricular assist device (LVAD)			
<b>Manufacturer name</b>			
Abiomed			
<b>Supplier name</b>			
Gada Italy			
<b>RDM</b>	<b>REF</b>		
129812	0550-0002		
<b>Type</b>	<b>CE marking (date)</b>	<b>Risk class</b>	<b>FDA approval</b>
2	First publication date in Ministry of Health database: 29/10/2020	III	Yes
<b>CND</b>			
Not applicable (assembled)			
<b>Clinical problem and rationale for the request</b>			
The Impella 5.5 has been discontinued and is replaced by the Impella 5.5 with Smart Assist at a price that has more than doubled (16,500 to 35,000 euros excluding VAT, respectively). In the new version, an extension has been added to the side arm of the bleeding tube of the Impella device.			
<b>Indications for use</b>			
<p>The Impella 5.5 Heart Pump with SmartAssist is an intracardiac pump for support of the left ventricle. It is intended for clinical use in cardiology and cardiac surgery for up to 30 days for the following non-exclusive indications</p> <ul style="list-style-type: none"> <li>-Impella 5.5 with SmartAssist heart pump is a cardiovascular support system for patients with reduced left ventricular function, including post cardiectomy, low output syndrome, cardiogenic shock after acute myocardial infarction.</li> <li>-It is also possible to use the Impella 5.5 heart pump with SmartAssist as a cardiovascular support system during beating-heart coronary bypass surgery, particularly in patients with a limited preoperative ejection fraction and a high risk of postoperative low output syndrome.</li> </ul>			
<b>Target patient</b>			
See Directions for Use.			

<b>General technology data</b>
<b>Main competitors</b>
Not applicable
<b>Technological details</b>
<b>Description</b>
A minimally invasive heart pump that provides peak flows of up to 5.5 l/min antegrade flow, enabling cardiac recovery. Impella 5.5 with SmartAssist can be inserted through the valve into the left ventricle with access via the axillary artery. It directly unloads the left ventricle by reducing ventricular work for up to 30 days for indications including support during cardiogenic shock. Provides peak flow of up to 5.5 l/min for systemic perfusion. Supports coronary perfusion. Allows repositioning in ICU without imaging (only for ventricle-migrated pumps). 70 cm catheter length facilitates patient mobilisation.
<b>Elements of Innovation</b>
Compared to the discontinued Impella 5.0 version, Impella 5.5 with Smart Assist features an extension to the side arm of the bleeding tube of the Impella device. It remains to be seen whether this change impacts on the efficacy and safety profile (see section 'Clinical studies').
<b>Clinical and economic evidence</b>
<b>Clinical studies</b>
The literature search selects a single clinical trial [1] whose objective was to evaluate whether, compared to standard practice, the addition of an extension to the side arm of the Impella device purge tube reduces the incidence of device failure and positively affects the clinical outcome of adult patients receiving Impella support. For this purpose, a retrospective review of the medical records of ICU patients in a tertiary care centre from August 2018 to August 2022 was performed to assess differences in patient outcomes related to the use of the Impella device before and after the implementation of the purge tube side arm extension. Twenty patients were included in the study of which 7 seven did not have the purge tube side arm extension, while 13 patients had the extension. The two groups showed no significant differences in clinical outcomes. In addition, there were no cases of failure and no cases of rupture of the purge cassette with the addition of the extension tube. There was a complete reduction in the incidence of purge cassette cracks, which could reduce infections. The authors report that long-term prospective studies are needed to confirm the results.
Pubmed search ( <a href="https://pubmed.ncbi.nlm.nih.gov/">https://pubmed.ncbi.nlm.nih.gov/</a> ): date: 31/05/2024; keyword: 'impella 5.5 smart assist'.
<b>Ongoing clinical trials</b>
Three studies are registered on ClinicalTrials.gov ( <a href="https://www.clinicaltrials.gov/">https://www.clinicaltrials.gov/</a> , search date: 05/06/2024, keyword: Impella 5.5):
<ol style="list-style-type: none"> <li>5. IMPACT EU Post-Market Clinical Follow-Up Study</li> <li>6. ClinicalTrials.gov ID NCT05756751. European, prospective, multicentre, single-arm, post-marketing clinical follow-up study. Up to 123 patients followed for up to 1 year will be enrolled in the study. The objective of the study is to collect data on the safety and efficacy of the Impella 5.5 system in alleviating post-cardiotomy heart failure and improving functional status and quality of life in high-risk patients undergoing cardiac surgery with severe LV dysfunction in a post-market setting.</li> <li>7. IMPella-Protected cArdiaC Surgery Trial (IMPACT). ClinicalTrials.gov ID NCT05529654. The aim of this study is to demonstrate the safety and efficacy of using Impella 5.5 in high-risk cardiac surgery patients, with the overall objective of evaluating whether the use of Impella 5.5 with SmartAssist (Impella 5.5) in the peri-operative phase improves early haemodynamics, end-organ function and clinical outcomes in patients with severely reduced LV function undergoing cardiac surgery. The study plans to enrol 100 patients.</li> </ol>

8. Impella-supported Off-pump Coronary Artery Bypass Grafting in High Risk Revascularizations: A Single Center Prospective Observational Study. ClinicalTrials.gov ID NCT04339452. The study plans to enrol 30 patients.

**Guidelines**

According to the ESC 2023 guidelines for the treatment of acute coronary syndromes, the role of mechanical circulatory support devices (MCS) (extracorporeal veno-arterial membrane oxygenation [VA-ECMO], micro-axial pumps) in the context of acute myocardial infarction is not well established and large-scale randomised studies are needed [2].

**Cost-effectiveness analysis**

-

**HTA Report**

-

**Expected benefits**

-

**Price and cost of therapy per patient**

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)	Annual requirements (No. pieces)
Impella 5.5 Smart Assist (Abiomed)	35,000 + 4% VAT	35.000	-

**Price and cost of therapy per patient with existing treatment alternatives**

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)
Impella 5.0* (Abiomed)	16,500 + 4% VAT	16,500 + 4% VAT

\*no longer on the market.

**Economic and organisational impact**

The use of Impella 5.5 with Smart Assist entails an increase in expenditure as this new version has a much higher price than the one no longer on the market.

**Reimbursability information**

Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM Code )	DRG code	Tariff (euro)
785.51 (Cardiogenic shock)	00.66 (Percutaneous Transluminal Coronary Angioplasty, PTCA, or coronary atherectomy)	127 (Heart failure and shock)	2.871

**Assessment of innovativeness (according to Regional Resolution N° 737/2022, [link](#) and N° 1244/2022 [link](#) )**

Innovative device (Y/N)	N
If yes, please indicate which Criteria 1, 2 and 3 are met:	

Summary data		
Request number	Date of request	Applicant
324	30/05/2024	Cardiac Surgery - AOUC
Technology		
IMPELLA CATHETER KIT 5.5 WITH SMART ASSIST (minimally invasive LVAD for the treatment of cardiogenic shock)		
Conclusions		
<p>Impella 5.5 with Smart Assist is a minimally invasive left ventricular assist device for the treatment of ongoing cardiogenic shock for up to 14 days. Impella aims to reduce the ventricular workload and provide the circulatory support necessary for myocardial recovery.</p> <p>Impella 5.5 with Smart Assist replaces the Impella 5.0 version, which has been withdrawn from the market. Impella 5.5 with Smart Assist was preliminarily evaluated in a retrospective study of 20 patients, the results of which are to be confirmed in larger studies [1].</p> <p>The price of this new version is considerably higher than that of the previous version (35,000 versus 16,500 euros excluding VAT, respectively).</p> <p>The absence of evidence obviously raises questions of both efficacy and safety. Moreover, in the absence of clinical data, the price increase does not seem justified. Based on these considerations, the CO defers to the decision of the Regional Commission.</p>		
Date of report		
31/05/2024		
Author(s) of the report		
Sabrina Trippoli		
Referring pharmacist		
Monica Vaiani		

## BIBLIOGRAPHY

1. Satashia P, White A, Isha S, Hanson A, Jenkins A, Blasavage J, Matos N, Tomlinson A, Zhang S, Quinones Q, Waldron N, Bhattacharyya A, Kiley S, Guru P, Chaudhury S, Shapiro A, Moreno Franco P, Sanghavi DK. Adding an extension piece to the end of the purge side arm of the Impella device can prevent the incidence of the cassette breaking and decrease the Impella device failure rate: Impact of practice change on patient outcome. *Perfusion*. 2023 Dec 7:2676591231220305. doi: 10.1177/02676591231220305. Epub ahead of print. PMID: 38060246.
2. Linee guida ESC 2023 per il trattamento delle sindromi coronariche acute. Indirizzo web: [https://www.giornaledicardiologia.it/r.php?v=4191&a=41785&l=354475&f=allegati/04191\\_2024\\_02/fulltext/04\\_LG\\_ESC\\_ACS%20\(e1-e112\).pdf](https://www.giornaledicardiologia.it/r.php?v=4191&a=41785&l=354475&f=allegati/04191_2024_02/fulltext/04_LG_ESC_ACS%20(e1-e112).pdf).