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2. Otto CM, Nishimura RA, Bonow RO, Carabello BA, Erwin JP 3rd, Gentile F, Jneid H, Krieger EV, Mack M, McLeod C, O'Gara PT, Rigolin VH, Sundt TM 3rd, Thompson A, Toly C; ACC/AHA Joint Committee Members; O'Gara PT, Beckman JA, Levine GN, Al-Khatib SM, Armbruster A, Birtcher KK, Cigarroa J, Deswal A, Dixon DL, Fleisher LA, de Las Fuentes L, Gentile F, Goldberger ZD, Gorenek B, Haynes N, Hernandez AF, Hlatky MA, Joglar JA, Jones WS, Marine JE, Mark D, Palaniappan L, Piano MR, Spatz ES, Tamis-Holland J, Wijeyesundera DN, Woo YJ. 2020 ACC/AHA guideline for the management of patients with valvular heart disease: A report of the American College of Cardiology/American Heart Association Joint Committee on Clinical Practice Guidelines. *J Thorac Cardiovasc Surg.* 2021 Aug;162(2):e183-e353. doi: 10.1016/j.jtcvs.2021.04.002. Epub 2021 May 8. PMID: 33972115.
3. Yu J, Qiao E, Wang W. Mechanical or biologic prostheses for mitral valve replacement: A systematic review and meta-analysis. *Clin Cardiol.* 2022 Jul;45(7):701-716. doi: 10.1002/clc.23854. Epub 2022 Jun 5. PMID: 35665516; PMCID: PMC9286334.
4. Stocco F, Fabozzo A, Bagozzi L, Cavalli C, Tarzia V, D'Onofrio A, Lorenzoni G, Chiminazzo V, Gregori D, Gerosa G. Biological versus mechanical aortic valve replacement in non-elderly patients: a single-centre analysis of clinical outcomes and quality of life. *Interact Cardiovasc Thorac Surg.* 2021 Apr 19;32(4):515-521. doi: 10.1093/icvts/ivaa306. PMID: 33313767; PMCID: PMC8906774.
5. Chikwe J, Lee T, Itagaki S, Adams DH, Egorova NN. Long-Term Outcomes After Off-Pump Versus On-Pump Coronary Artery Bypass Grafting by Experienced Surgeons. *J Am Coll Cardiol.* 2018 Sep 25;72(13):1478-1486. doi: 10.1016/j.jacc.2018.07.029. PMID: 30236310.
6. Vahanian A, Beyersdorf F, Praz F, Milojevic M, Baldus S, Bauersachs J, Capodanno D, Conradi L, De Bonis M, De Paulis R, Delgado V, Freemantle N, Gilard M, Haugaa KH, Jeppsson A, Jüni P, Pierard L, Prendergast BD, Sádaba JR, Tribouilloy C, Wojakowski W; ESC/EACTS Scientific Document Group. 2021 ESC/EACTS Guidelines for the management of valvular heart disease. *Eur Heart J.* 2022 Feb 12;43(7):561-632. doi: 10.1093/eurheartj/ehab395. Erratum in: *Eur Heart J.* 2022 Feb 18;: PMID: 34453165.
7. Shang H, Claessens SM, Tian B, Wright GA. Aldehyde reduction in a novel pericardial tissue reduces calcification using rabbit intramuscular model. *J Mater Sci Mater Med.* 2017 Jan;28(1):16. doi: 10.1007/s10856-016-5829-8. Epub 2016 Dec 20. PMID: 28000112; PMCID: PMC5174141.
8. Heimansohn DA, Baker C, Rodriguez E, Takayama H, Dagenais F, Talton DS, Mumtaz MA, Pibarot P, Puskas JD; COMMENCE Trial Investigators. Mid-term outcomes of the COMMENCE trial investigating mitral valve replacement using a bioprosthesis with a novel tissue. *JTCVS Open.* 2023 Jun 2;15:151-163. doi: 10.1016/j.xjon.2023.05.008. Erratum in: *JTCVS Open.* 2023 Dec 27;17:184. PMID: 37808026; PMCID: PMC10556809.
9. Imada R, Nagatomi S, Kinjo T, Oono H, Tateishi N, Kawazu Y. [Transprosthetic Cuff Leakage of the MITRIS RESILIA Mitral Valve During Mitral Valve Replacement: Report of a Case]. *Kyobu Geka.* 2022 Dec;75(13):1112-1116. Japanese. PMID: 36539229.
10. Heye T, Reemtsen B, Greiten L. The MITRIS RESILIA mitral valve is a safe and effective option for mitral valve replacement in young patients requiring mitral valve replacement. *Cardiol Young.* 2022 Sep 28:1-3. doi: 10.1017/S1047951122003043. Epub ahead of print. PMID: 36169001.
11. Ushijima T, Kimura S, Shiose A. Implantability of the MITRIS RESILIA Mitral Valve. *Ann Thorac Surg.* 2022 Apr;113(4):e295-e297. doi: 10.1016/j.athoracsur.2021.05.065. Epub 2021 Jun 17. PMID: 34147490.
12. Bavaria JE, Griffith B, Heimansohn DA, Rozanski J, Johnston DR, Bartus K, Girardi LN, Beaver T, Takayama H, Mumtaz MA, Rosengart TK, Starnes V, Timek TA, Boateng P, Ryan W, Cornwell LD, Blackstone EH, Borger MA, Pibarot P, Thourani VH, Svensson LG, Puskas JD; COMMENCE Trial Investigators. Five-year Outcomes of the COMMENCE Trial Investigating Aortic Valve Replacement With RESILIA Tissue. *Ann Thorac Surg.* 2023 Jun;115(6):1429-1436. doi: 10.1016/j.athoracsur.2021.12.058. Epub 2022 Jan 20. PMID: 35065065.
13. Beute TJ, Goehler M, Parker J, Boeve T, Heiser J, Murphy E, Timek T, Willekes CL. Long-Term Outcomes of Mosaic Versus Perimount Mitral Replacements: 17-Year Follow-Up of 940 Implants. *Ann Thorac Surg.* 2020 Aug;110(2):508-515. doi: 10.1016/j.athoracsur.2019.10.075. Epub 2019 Dec 20. PMID: 31866480.
14. Loor G, Schuster A, Cruz V, Rafael A, Stewart WJ, Diaz J, McCurry K. The Carpentier-Edwards Perimount Magna mitral valve bioprosthesis: intermediate-term efficacy and durability. *J Cardiothorac Surg.* 2016 Jan 27;11:20. doi: 10.1186/s13019-016-0412-4. PMID: 26818795; PMCID: PMC4728780.

15. Chiariello GA, Beraud AS, Vahdat O, Van Rothem J, Garcia O, Soula P, Berthoumieu P, Abouliatim I. Late results after mitral valve replacement with Mosaic bioprosthesis in patients aged 65 years or younger. *Interact Cardiovasc Thorac Surg.* 2021 Jul 26;33(2):181-187. doi: 10.1093/icvts/ivab066. PMID: 33693682; PMCID: PMC8923380.
16. Nakazato T, Hata H, Toda K, Miyagawa S, Yoshikawa Y, Saito S, Domae K, Ueno T, Kuratani T, Sawa Y. Midterm Clinical Outcomes of the St Jude Medical Epic Porcine Bioprosthesis in the Mitral Position. *Circ J.* 2018 Dec 25;83(1):110-116. doi: 10.1253/circj.CJ-18-0483. Epub 2018 Oct 25. PMID: 30369546.
17. Zwischenberger BA, Gaca JG, Haney JC, Carr K, Glower DD. Late Outcomes of Porcine and Pericardial Bioprostheses after Mitral Valve Replacement in 1162 Patients. *Ann Thorac Surg.* 2024 Apr 15:S0003-4975(24)00286-8. doi: 10.1016/j.athoracsur.2024.04.002. Epub ahead of print. PMID: 38631662.
18. Koulouroudias M, Di Mauro M, Chiariello G, Meani P, Lorusso R. Long-Term Outcomes of Bioprosthetic Valves in the Mitral Position: A Pooled Meta-Analysis of Reconstructed Time-to-Event Individual Patient Data. *Am J Cardiol.* 2024 Apr 16;221:64-73. doi: 10.1016/j.amjcard.2024.04.008. Epub ahead of print. PMID: 38636624.

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RAPID HTA REPORT		
Request No.	Date of request	Applicant
322	02/05/2024	Cardiac Surgery - AOUC
Type of report		
New report		X
Updating a previous report		
If updating, please state the reason:		

General technology data					
<b>Trade name</b>					
Mitris Resilia					
<b>Generic name</b>					
MITRIS RESILIA Mitral Valve					
<b>Manufacturer name</b>					
Edwards Lifesciences					
<b>Supplier name</b>					
Edwards Lifesciences					
<b>RDM</b>	<b>REF</b>				
2491202	11400M25				
2491253	11400M27				
2491258	11400M29				
2491262	11400M31				
2491267	11400M33				
<b>Type</b>	<b>CE marking (date)</b>	<b>Risk class</b>	<b>FDA approval</b>		
1	12/10/2023	III	Yes		
<b>CND</b>					
P070301202 - BIOLOGICAL MITRAL VALVES WITH SURGICAL IMPLANT SUPPORT, NON-VALVULAR TISSUE OF ANIMAL ORIGIN					
<b>Clinical problem and rationale for the request</b>					
<p>Mitral valve disease is one of the most common valve diseases with an increasing prevalence in developed countries due to an ageing population [1]. According to heart valve disease guidelines, mitral valve replacement (MVR) may be considered in patients with mitral regurgitation and mitral stenosis when a durable repair is not feasible [2]. The optimal prosthetic valve is still a matter of debate, particularly for patients younger than 65 years of age [3]. Patients can be implanted with a bioprosthetic valve or a mechanical valve. Mechanical valves may offer longer durability, but require lifelong anticoagulant therapy with a potential risk of haemorrhage [4]. In contrast, bioprosthetic valves have been associated with an increased risk of structural valve deterioration especially in younger patients and reduced survival in patients younger than 50 years [5]. The European Society of Cardiology suggests intervention with mechanical valves for patients younger than 65 years and with bioprostheses for patients older than 70 years [6]. Among these ages, recent consensus guidelines recommend the choice of prosthetic valves based on patient preference and expert opinion, and suggest that it should be a shared decision-making process.</p> <p>RESILIA tissue, designed for use in heart valve replacement, is a promising option that can improve the durability of bioprosthetic valves. This tissue has shown reduced calcification in preclinical and animal studies compared to other tissues [7].</p>					

General technology data
<b>Indications for use</b>
The MITRIS RESILIA mitral valve is intended for use as a heart valve replacement. It is indicated for replacement of the native or prosthetic mitral valve
<b>Target patient</b>
Young subjects under 75 years of age with mitral steno-insufficiency who are not candidates for valve repair in the setting of rheumatic disease and/or new valve degeneration following radiation therapy. Severe calcifications of the valvular apparatus leading to complex steno-deficiency conditions are observed in these cases.
<b>Main competitors</b>
Magna Mitral Ease Valve (mitral tissue biological valve of animal origin - bovine pericardium - with support awarded to Lot 3 of Determination No. 1402/2022 effective 1.1.2023)

Technological details
<b>Description</b>
Mitris Resilia is a three-flap tissue valve made of animal tissue (RESILIA bovine pericardium) and fitted with a nitinol stent; a cobalt-chrome strip and a polyester strip surround the base of the valve under the wire mesh frame providing structural support for the orifice. The two metal alloys are visible on fluoroscopy for identification of the inflow and outflow margins. A polytetrafluoroethylene (PTFE) suture ring is attached to the wire mesh frame to facilitate tissue integration and valve orientation at the time of placement.
<b>Elements of Innovation</b>
The Resilia fabric of which Mitris valves are made incorporates an anti-calcification process that blocks residual aldehyde groups that bind to calcium. The technology includes tissue preservation in glycerol that eliminates tissue exposure to unbound residual aldehyde groups commonly found in glutaraldehyde solutions. Resilia tissue also allows valves to be stored dry rather than in the storage fluid, simplifying and speeding up work in the operating room. Recent data from the COMMENCE study [8], see section "Clinical studies", show that Resilia tissue, 5 years after implantation, results in low rates of structural valve deterioration, often caused by calcium build-up, which can affect the long-term durability of bioprosthetic valves. The Mitris valve shows good results in terms of safety and effectiveness. Currently implanted bioprosthetic valves in the mitral position have an average durability of about 7 years, and this is reduced in patients under 65 years of age.

Clinical and economic evidence
<b>Clinical studies</b>
The literature search selected one observational study (Commence trial limited to patients undergoing mitral valve replacement [8]), and three case reports [9-11].
The COMMENCE (FAD Investigational Device Exemption) trial was planned to evaluate the performance of RESILIA in aortic and mitral valve replacement [11]. The study described here reports the mid-term results of the COMMENCE trial of the mitral cohort [8].
Between January 2013 and February 2016, 83 patients were enrolled at 17 centres in North America, with 82 patients successfully implanted with the study valve. 54 patients had data available at 5-year follow-up. The patients were predominantly elderly and had concomitant cardiovascular disease. The average age of the patients at the time of surgery was 70 years, with 52.4 % of patients aged 70 years or older.
The diagnoses were as follows: 48.8% pure insufficiency (n = 40/82), 30.5% stenosis with insufficiency (n = 25/82), 11.0% stenosis (n = 9/82), 7.3% prosthetic valve dysfunction (n = 6/82) and 2.4% other diagnosis (n = 2/82). The majority of patients were taking antiplatelet therapy (51.2%) and 46.3% of patients were taking anticoagulant therapy at baseline. The all-cause mortality within 30 days after surgery was 1.2%. There were 2 procedure-related

**Clinical and economic evidence**

early ischaemic strokes (2.4%) and 1 early haemorrhagic event (1.2%) due to an oesophageal laceration, which was also considered procedure-related. No valve explants, valve thrombosis, endocarditis, haemolysis or reinterventions were observed in the 30 days following the procedure. At 5-year follow-up, freedom from all-cause mortality, structural valve deterioration and re-intervention was 79.9%, 98.7% and 97.1%, respectively. Haemodynamic measures of valve function remained stable during the 5-year follow-up period; valve leakage was observed infrequently and mostly clinically insignificant/mild.

At 5-year follow-up, the COMMENCE Mitral study reported a good safety profile of the Mitral Mitris valve with RESILIA tissue. These results are consistent with the mid-term safety and efficacy results observed with RESILIA tissue in the aortic position [12].

A comparative analysis of 940 patients undergoing surgery with the Medtronic Mosaic or PERIMOUNT mitral valves reported 5-year freedom from structural valve deterioration rates of 91.0% and 90.3%, respectively [13]. Furthermore, publications on the PERIMOUNT Magna, Medtronic Mosaic and Abbott Epic bioprosthetic valves reported 5-year freedom from structural valve deterioration rates of 90%, 94% and 93%, respectively [14-16]. The freedom from structural deterioration rate observed in the COMMENCE trial is in line with that of the other already published mid-term studies [14-16].

A limitation of the study is the small size of the patient cohort. Furthermore, as several studies show a significant increase in valve deterioration from 5 years onwards [13-16], long-term follow-up is necessary to understand whether these results are maintained beyond 5 years. For this reason, the results of the extended 10-year follow-up cohort of the COMMENCE study are important, as are the results of the ongoing MOMENTIS prospective trial (ClinicalTrials.gov, NCT05526560), which collects data on the acute and long-term safety and performance of the MITRIS RESILIA mitral valve, with results reported 10 years after implantation.

Pubmed search (<https://pubmed.ncbi.nlm.nih.gov/>): keyword 'mitris resilia'; date 15/05/2024.

**Ongoing clinical trials**

The trial 'Real-world Clinical Outcomes of the MITRIS RESILIA Mitral Valve (MOMENTIS)' is registered on ClinicalTrials.

(NCT05526560). MOMENTIS is a prospective, observational, single-arm, multicentre study designed to collect clinical outcomes in up to 500 subjects who received the MITRIS RESILIA mitral valve, model 11400M for native or prosthetic mitral valve replacement with or without concomitant procedures. The study is scheduled to end in December 2035. The primary end-points are freedom from death or valve reintervention at 12 months after surgery and haemodynamic performance as measured by peak and mean gradient by echocardiography at 12 months after surgery.

Search on ClinicalTrials.gov (<https://clinicaltrials.gov/>): keywords "mitris resilia", date 15/05/2024.

**Guidelines**

There are no guidelines as to which type of biological valve (porcine or bovine origin) should be used. Two very recent publications found that there is no difference in survival time and structural deterioration between bovine and porcine pericardium valves in the long term. Structural deterioration is associated with the younger age of the patient, but not with the type of valve [17-18].

**Cost-effectiveness analysis**

No cost-effectiveness analysis available

**HTA Report**

No HTA reports available

**Expected benefits**

**Clinical and economic evidence**

Increased valve life.

**Price and cost of therapy per patient**

Product (Manufacturer)	Unit price* (euro)	Therapy cost per patient (euro)	Annual requirements (No. pieces)
Mitris Resilia (Edwards Lifesciences)	3.850	3.850	30

\*VAT excluded at 4%.

**Price and cost of therapy per patient with existing treatment alternatives**

Product (Manufacturer)	Unit price* (euro)	Therapy cost per patient (euro)
Mitral Mosaic <sup>^</sup> (Medtronic)	1.550	1.550
EPIC SUPRA <sup>^</sup> (Abbott)	1.400	1.400
MAGNA MITRAL EASE <sup>^^</sup> (Edwards Lifesciences)	2.500	2.500

\*VAT excluded at 4%.

<sup>^</sup> biologic mitral heart valves from porcine valve tissue with support awarded Lot 1 of Determination No. 1402/2022 effective 1.1.2023. In 2023 Mosaic Mitral was purchased by all three AOU and the Monasterio Foundation. Epic Supra was purchased by the AOUP and the Monasterio Foundation (source: BO, data as at 12/05/2024).

<sup>^^</sup>Biological mitral heart valve from animal tissue (bovine pericardium) with support awarded in Lot 3 of Determination No. 1402/2022 effective 1.1.2023. This valve was purchased by all three AOU (source: Qlik, data as at 12/05/2024).

**Economic and organisational impact**

The use of this valve entails an increase in expenditure as it has a higher price than the comparator, Magna Mitral Ease valve, in use.

**Reimbursability information**

Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM Code )	DRG code	Tariff (euro)
394.0 - 394.2 - 396.0 - 394.1 - 396.1	35.23	-	-

**Assessment of innovativeness (according to Regional Resolution N° 737/2022, [link](#) and N° 1244/2022 [link](#) )**

Innovative device (Y/N)	N
If yes, please indicate which Criteria 1, 2 and 3 are met:	

**Summary data**

Request number	Date of request	Applicant
322	02/05/2024	Cardiac Surgery - AOUC
Technology		
MITRIS RESILIA mitral valve (biological mitral valve from bovine pericardium for the treatment of mitral stenosis in patients who are not candidates for valve repair)		
Conclusions		

**Summary data**

The Mitris Resilia Valve from Edwards Lifesciences recently received CE marking for replacement of the native or prosthetic mitral valve. It is made from bovine pericardial tissue (Resilia) which, by reducing calcium build-up on the valve, should preserve and protect the tissue, allowing the valve to potentially last longer than traditional bioprosthetic valves whose average durability is about 7 years, which is reduced in patients under 65 years of age [6].

Recent data from the COMMENCE study [8] show that Resilia tissue results in low rates of structural valve deterioration at a 5-year follow-up after implantation. These results are in line with those observed with other animal-derived valves [13-16]. However, as the literature shows a significant increase in deterioration of tissue valves from 5 years onwards [13-16], long-term follow-up is necessary to understand whether the results of the COMMENCE trial are maintained beyond 5 years. The outcomes of the extended 10-year follow-up cohort of the COMMENCE trial and the outcomes of the ongoing prospective MOMENTIS trial (ClinicalTrials.gov. NCT05526560), with results at 10 years after implantation, may define this.

Valves of animal origin for mitral valve replacement (Mosaic Mitral, Epic Supra and Magna Mitral Ease) are already available in Tuscany. To date, it is not known whether Mitris Resilia has any advantages, especially in terms of structural deterioration of the valve, in comparison with the valves already in use since, as mentioned above, long-term data are lacking. Nevertheless, Mitris Resilia is priced higher than the comparators.

Mitris Resilia did not obtain recognition of innovativeness on the basis of the criteria defined by the Region of Tuscany (Resolution No. 737/2022).

Based on these considerations, the CO expresses an unfavourable opinion on the purchase.

**Date of report**

20/05/2024

**Author(s) of the report**

Sabrina Trippoli

**Referring pharmacist**

Monica Vaiani

**BIBLIOGRAPHY**

19. Tsao CW, Aday AW, Almarzooq ZI, Anderson CAM, Arora P, Avery CL, Baker-Smith CM, Beaton AZ, Boehme AK, Buxton AE, Commodore-Mensah Y, Elkind MSV, Evenson KR, Eze-Nliam C, Fugar S, Generoso G, Heard DG, Hiremath S, Ho JE, Kalani R, Kazi DS, Ko D, Levine DA, Liu J, Ma J, Magnani JW, Michos ED, Mussolino ME, Navaneethan SD, Parikh NI, Poudel R, Rezk-Hanna M, Roth GA, Shah NS, St-Onge MP, Thacker EL, Virani SS, Voeks JH, Wang NY, Wong ND, Wong SS, Yaffe K, Martin SS; American Heart Association Council on Epidemiology and Prevention Statistics Committee and Stroke Statistics Subcommittee. Heart Disease and Stroke Statistics-2023 Update: A Report From the American Heart Association. *Circulation*. 2023 Feb 21;147(8):e93-e621. doi: 10.1161/CIR.0000000000001123. Epub 2023 Jan 25. Erratum in: *Circulation*. 2023 Feb 21;147(8):e622. Erratum in: *Circulation*. 2023 Jul 25;148(4):e4. PMID: 36695182.
20. Otto CM, Nishimura RA, Bonow RO, Carabello BA, Erwin JP 3rd, Gentile F, Jneid H, Krieger EV, Mack M, McLeod C, O'Gara PT, Rigolin VH, Sundt TM 3rd, Thompson A, Toly C; ACC/AHA Joint Committee Members; O'Gara PT, Beckman JA, Levine GN, Al-Khatib SM, Armbruster A, Birtcher KK, Cigarroa J, Deswal A, Dixon DL, Fleisher LA, de Las Fuentes L, Gentile F, Goldberger ZD, Gorenek B, Haynes N, Hernandez AF, Hlatky MA, Joglar JA, Jones WS, Marine JE, Mark D, Palaniappan L, Piano MR, Spatz ES, Tamis-Holland J, Wijeyesundera DN, Woo YJ. 2020 ACC/AHA guideline for the management of patients with valvular heart disease: A report of the American College of Cardiology/American Heart Association Joint Committee on Clinical Practice Guidelines. *J Thorac Cardiovasc Surg*. 2021 Aug;162(2):e183-e353. doi: 10.1016/j.jtcvs.2021.04.002. Epub 2021 May 8. PMID: 33972115.

21. Yu J, Qiao E, Wang W. Mechanical or biologic prostheses for mitral valve replacement: A systematic review and meta-analysis. *Clin Cardiol.* 2022 Jul;45(7):701-716. doi: 10.1002/clc.23854. Epub 2022 Jun 5. PMID: 35665516; PMCID: PMC9286334.
22. Stocco F, Fabozzo A, Bagozzi L, Cavalli C, Tarzia V, D'Onofrio A, Lorenzoni G, Chiminazzo V, Gregori D, Gerosa G. Biological versus mechanical aortic valve replacement in non-elderly patients: a single-centre analysis of clinical outcomes and quality of life. *Interact Cardiovasc Thorac Surg.* 2021 Apr 19;32(4):515-521. doi: 10.1093/icvts/ivaa306. PMID: 33313767; PMCID: PMC8906774.
23. Chikwe J, Lee T, Itagaki S, Adams DH, Egorova NN. Long-Term Outcomes After Off-Pump Versus On-Pump Coronary Artery Bypass Grafting by Experienced Surgeons. *J Am Coll Cardiol.* 2018 Sep 25;72(13):1478-1486. doi: 10.1016/j.jacc.2018.07.029. PMID: 30236310.
24. Vahanian A, Beyersdorf F, Praz F, Milojevic M, Baldus S, Bauersachs J, Capodanno D, Conradi L, De Bonis M, De Paulis R, Delgado V, Freemantle N, Gilard M, Haugaa KH, Jeppsson A, Jüni P, Pierard L, Prendergast BD, Sádaba JR, Tribouilloy C, Wojakowski W; ESC/EACTS Scientific Document Group. 2021 ESC/EACTS Guidelines for the management of valvular heart disease. *Eur Heart J.* 2022 Feb 12;43(7):561-632. doi: 10.1093/eurheartj/ehab395. Erratum in: *Eur Heart J.* 2022 Feb 18;: PMID: 34453165.
25. Shang H, Claessens SM, Tian B, Wright GA. Aldehyde reduction in a novel pericardial tissue reduces calcification using rabbit intramuscular model. *J Mater Sci Mater Med.* 2017 Jan;28(1):16. doi: 10.1007/s10856-016-5829-8. Epub 2016 Dec 20. PMID: 28000112; PMCID: PMC5174141.
26. Heimansohn DA, Baker C, Rodriguez E, Takayama H, Dagenais F, Talton DS, Mumtaz MA, Pibarot P, Puskas JD; COMMENCE Trial Investigators. Mid-term outcomes of the COMMENCE trial investigating mitral valve replacement using a bioprosthesis with a novel tissue. *JTCVS Open.* 2023 Jun 2;15:151-163. doi: 10.1016/j.jxon.2023.05.008. Erratum in: *JTCVS Open.* 2023 Dec 27;17:184. PMID: 37808026; PMCID: PMC10556809.
27. Imada R, Nagatomi S, Kinjo T, Oono H, Tateishi N, Kawazu Y. [Transprosthetic Cuff Leakage of the MITRIS RESILIA Mitral Valve During Mitral Valve Replacement: Report of a Case]. *Kyobu Geka.* 2022 Dec;75(13):1112-1116. Japanese. PMID: 36539229.
28. Heye T, Reemtsen B, Greiten L. The MITRIS RESILIA mitral valve is a safe and effective option for mitral valve replacement in young patients requiring mitral valve replacement. *Cardiol Young.* 2022 Sep 28:1-3. doi: 10.1017/S1047951122003043. Epub ahead of print. PMID: 36169001.
29. Ushijima T, Kimura S, Shiose A. Implantability of the MITRIS RESILIA Mitral Valve. *Ann Thorac Surg.* 2022 Apr;113(4):e295-e297. doi: 10.1016/j.athoracsur.2021.05.065. Epub 2021 Jun 17. PMID: 34147490.
30. Bavaria JE, Griffith B, Heimansohn DA, Rozanski J, Johnston DR, Bartus K, Girardi LN, Beaver T, Takayama H, Mumtaz MA, Rosengart TK, Starnes V, Timek TA, Boateng P, Ryan W, Cornwell LD, Blackstone EH, Borger MA, Pibarot P, Thourani VH, Svensson LG, Puskas JD; COMMENCE Trial Investigators. Five-year Outcomes of the COMMENCE Trial Investigating Aortic Valve Replacement With RESILIA Tissue. *Ann Thorac Surg.* 2023 Jun;115(6):1429-1436. doi: 10.1016/j.athoracsur.2021.12.058. Epub 2022 Jan 20. PMID: 35065065.
31. Beute TJ, Goehler M, Parker J, Boeve T, Heiser J, Murphy E, Timek T, Willekes CL. Long-Term Outcomes of Mosaic Versus Perimount Mitral Replacements: 17-Year Follow-Up of 940 Implants. *Ann Thorac Surg.* 2020 Aug;110(2):508-515. doi: 10.1016/j.athoracsur.2019.10.075. Epub 2019 Dec 20. PMID: 31866480.
32. Loor G, Schuster A, Cruz V, Rafael A, Stewart WJ, Diaz J, McCurry K. The Carpentier-Edwards Perimount Magna mitral valve bioprosthesis: intermediate-term efficacy and durability. *J Cardiothorac Surg.* 2016 Jan 27;11:20. doi: 10.1186/s13019-016-0412-4. PMID: 26818795; PMCID: PMC4728780.
33. Chiariello GA, Beraud AS, Vahdat O, Van Rothem J, Garcia O, Soula P, Berthoumieu P, Abouliatim I. Late results after mitral valve replacement with Mosaic bioprosthesis in patients aged 65 years or younger. *Interact Cardiovasc Thorac Surg.* 2021 Jul 26;33(2):181-187. doi: 10.1093/icvts/ivab066. PMID: 33693682; PMCID: PMC8923380.
34. Nakazato T, Hata H, Toda K, Miyagawa S, Yoshikawa Y, Saito S, Domae K, Ueno T, Kuratani T, Sawa Y. Midterm Clinical Outcomes of the St Jude Medical Epic Porcine Bioprosthesis in the Mitral Position. *Circ J.* 2018 Dec 25;83(1):110-116. doi: 10.1253/circj.CJ-18-0483. Epub 2018 Oct 25. PMID: 30369546.
35. Zwischenberger BA, Gaca JG, Haney JC, Carr K, Glower DD. Late Outcomes of Porcine and Pericardial Bioprostheses after Mitral Valve Replacement in 1162 Patients. *Ann Thorac Surg.* 2024 Apr 15:S0003-4975(24)00286-8. doi: 10.1016/j.athoracsur.2024.04.002. Epub ahead of print. PMID: 38631662.

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***Centro operativo***

36. Koulouroudias M, Di Mauro M, Chiariello G, Meani P, Lorusso R. Long-Term Outcomes of Bioprosthetic Valves in the Mitral Position: A Pooled Meta-Analysis of Reconstructed Time-to-Event Individual Patient Data. *Am J Cardiol.* 2024 Apr 16;221:64-73. doi: 10.1016/j.amjcard.2024.04.008. Epub ahead of print. PMID: 38636624.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Reporting by the Operations Centre, Regional Decree No. 17610 of 7 September 2022.

For further information write to [centro.operativo.htart@regione.toscana.it](mailto:centro.operativo.htart@regione.toscana.it)