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Fabiola Del Santo

BIBLIOGRAFIA

1. Kusumoto FM, Schoenfeld MH, Wilkoff BL, Berul CI, Birgersdotter-Green UM, Carrillo R et al. 2017 HRS expert consensus statement on cardiovascular implantable electronic device lead management and extraction. *Heart Rhythm* 2017;14:e503–51.
2. Dai M, Cai C, Vaibhav V, Sohail MR, Hayes DL, Hodge DO et al. Trends of cardiovascular implantable electronic device infection in 3 decades: a population-based study. *JACC Clin Electrophysiol* 2019;5:1071–80.
3. Blomstrom-Lundqvist C, Ostrowska B. Prevention of cardiac implantable electronic device infections: guidelines and conventional prophylaxis. *Europace*. 2021 May 25;23(Suppl 4):iv11–9. doi: 10.1093/europace/euab071. Epub ahead of print. PMID: 34037227; PMCID: PMC8221047.
4. Shariff N, Akthar T, Razak E, Segerson N, Schwartzman D. Cardiovascular Implantable Electronic Device Infections: Risk Scoring and Role of Antibiotic Envelope in Prevention. *Recent Adv Cardiovasc Drug Discov*. 2015;10(1):70-6. doi: 10.2174/1574890110666151104120207. PMID: 26530457.
5. Clémenty N, Carion PL, Léotoing L, Lamarsalle L, Wilquin-Bequet F, Brown B, Verhees KJP, Fernandes J, Deharo JC. Infections and associated costs following cardiovascular implantable electronic device implantations: a nationwide cohort study. *Europace*. 2018 Dec 1;20(12):1974-1980. doi: 10.1093/europace/eux387. PMID: 29672690.
6. Borov S, Baldauf B, Henke J, Pavaci H, Perani A, Zrenner B, Dietl J, Mehilli J, Lau EW, Vonthein R, Bonnemeier H. Use of a taurolidine containing antimicrobial wash to reduce cardiac implantable electronic device infection. *Europace*. 2023 Oct 5;25(10):euad306. doi: 10.1093/europace/euad306. PMID: 37831737; PMCID: PMC10616572.
7. Giudice M, Catuzzo B, Berlier N, Lau EW, Bonnemeier H, Assadian O, Baldauf B, Borov S, Scacciarella P. Use of Taurolidine in a Patient With a Cardiac Implantable Electronic Device Protrusion. *JACC Case Rep*. 2023 Apr 19;14:101835. doi: 10.1016/j.jaccas.2023.101835. PMID: 37152697; PMCID: PMC10157148.
8. Giaccardi M, Baldauf B, Lau EW, Borov S, Bonnemeier H. Salvage of Cardiac Implantable Electronic Device Pocket Infection with Skin Erosion in Frail 92-Year-Old. *J Cardiovasc Dev Dis*. 2022 Mar 10;9(3):81. doi: 10.3390/jcdd9030081. PMID: 35323629; PMCID: PMC8955956.
9. British Heart Rhythm Society. Standards for implantation and follow-up of cardiac rhythm management devices in adults. January 2024 revision, <https://bhrc.com/wp-content/uploads/2024/01/BHRS-standards-January-2024-Implantation-and-Follow-Up-of-CRM-Devices-in-Adults.pdf>
10. Ludwig S, Theis C, Brown B, Witthohn A, Lux W, Goette A. Incidence and costs of cardiac device infections: retrospective analysis using German health claims data. *J Comp Eff Res*. 2018 May;7(5):483-492. doi: 10.2217/ce-2017-0080. Epub 2017 Nov 14. PMID: 29132224.
11. Love CJ. Taurolidine: could this be the 'silver bullet' against cardiac implantable electronic device infection? *Europace*. 2023 Nov 2;25(11):euad332. doi: 10.1093/europace/euad332. Erratum in: *Europace*. 2023 Dec 6;25(12): PMID: 37967367; PMCID: PMC10651068.

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RAPID HTA REPORT			
Request No.	Date of request	Applicant	
317	05/04/2024	UOC Cardiology Arezzo- UOS Interventional Arrhythmology	
Type of report			
New report			N
Update of a previous report (sheet 245 of 20/05/2022)			S
If update, please state reason: This sheet is an update of a previous evaluation (sheet 245 of 20/05/2022) on TauroPace. A trial has recently become available reporting a reduction in CIED-related infections with the use of TauroPace. TauroPace is mentioned in the British Heart Rhythm Society (BHRS) Guidelines as a possible treatment option to prevent CIED infections in patients at high risk of infection.			

General technology data					
Trade name					
TauroPace					
Generic name					
Antimicrobial solution for bacterial contamination CIED devices - 100 ml ampoules					
Manufacturer name					
TAUROPHARM GMBH					
Supplier name					
Medival srl					
RDM	REF				
1835797	TP-08-100				
Type	CE marking (date)	Risk class	FDA approval		
1	2174842CE02 (29/01/2019 date of CE marking and start of marketing).	IIb	No		
CND					
J010180 - PACE MAKER - ACCESSORIES					
Clinical problem and rationale for the request					
<p>Implantable cardiac electronic devices (CIEDs) are implanted in approximately 1.2-1.4 million patients per year worldwide. In recent years, an increase in the infection rate of CIEDs has been observed in relation to the expanding volume of procedures. This may be due to the increased number of co-morbidities presented by CIED patients, the increased complexity of implant systems, and the increase in life expectancy, which therefore leads to more revision procedures [1]. In a very long-term study [2], the cumulative odds of CIED infection reported were 6.2% at 15 years and 11.7% at 25 years, and 2.6%, 2.7% and 24.1% after the first, second and third procedures at 15 years. The hazard ratio (HR) for infection was 3.91 for pulse generator replacement and 3.08 for system upgrade. CIED infection generally requires complete explantation of all hardware. Lead extraction is a highly invasive procedure, associated with significant morbidity and potential mortality, and requires specialised training and equipment. Possible effective solutions against CIED infections currently are pre-operative antibiotics and antibiotic-eluting pockets, which have been tested in randomised clinical trials and are recommended in the current ESC guidelines [3].</p>					

General technology data

Indications for use

TauroPace is an extracorporeal disinfectant solution indicated for the eradication of environmental microbial contamination from the surface of CIEDs (implantable cardiac devices), prior to implantation in adult patients. Its application is therefore intended for any type of CIED, whether titanium or stainless steel, polyurethane or silicone, including its components (cables, leads, etc.).

Target patient

Patients implanted with CIEDs at high infectious risk. Patients at high infectious risk are those identified by risk scores specific to this population such as the Shariff et al. score [4] and the Padit score [5].

Main competitors

Antibiotic prophylaxis (mainly with cefazolin) is the standard for the prevention of CIED-related infections [3]. In addition to antibiotic prophylaxis, the TYRX device, a resorbable antibacterial wrap designed to contain a pacemaker or cardiac implantable defibrillator to reduce the risk of infection, is commercially available.

Technological details

Description

Tauropace is an antimicrobial solution for extracorporeal use with broad-spectrum antimicrobial properties for adult patients with implantable electronic cardiac devices (CIEDs). Tauropace consists of taurolidine in an aqueous solution (20 mg/ml); the active ingredient, an amine derived from the amino acid taurine, guarantees not only biocidal activity against fungi, Gram-negative and Gram-positive bacteria, preventing the formation of biofilm and subsequent infections, but also ensures long-lasting activity without the risk of resistance emergence. The surface of the CIEDs must be moistened with Tauropace to create an environment hostile to microbial adhesion and proliferation; this procedure, in order to ensure total coverage, must be performed not only on the generators but also on all components connected to them. CIEDs can be rubbed with gauze previously wetted with solution or directly immersed in the solution, poured in advance into a sterile container. The composition of the Tauropace solution does not create any interference with electronic parts and connection points

Elements of Innovation

Taurolidine with its lower cost and ease of application could be a cost-effective solution to CIED-related infections.

Clinical and economic evidence

Clinical studies

Borov's recent study [6] compared the efficacy of TauroPace with a control group treated with 3% hydrogen peroxide (H₂O₂) for the prevention of CIED infections. During CIED procedures, the device bag was irrigated and all accessible hardware was treated with H₂O₂ (n=532 patients), TauroPace or taurolidine in a galenic formulation (n=631 patients). The primary endpoint was major CIED infection within 3 months, while secondary endpoints included infections after 3 months, adverse events potentially related to antimicrobial solutions, CIED system, procedure and mortality until the end of follow-up. The TauroPace-treated group had more risk factors for infection than the H₂O group, but similar device- and procedure-specific risk factors. At 3 months, no CIED infection was observed in the TauroPace group (0/654, 0%), whereas 6 infections occurred in the H₂O₂ group (6/551, 1.1%), p = 0.0075. After the first 3 months, the TauroPace group had 3 infections and the H₂O₂ group had 9 infections (HR, 0.41; 95% CI: 0.10-1.56). Mortality was similar between the two groups (3.5% in the TauroPace group and 2.5% in the H₂O₂ group, p = 0.33). Non-infection-related adverse events were rarer in the TauroPace group (3.8%) than in the H₂O₂ group (6.0%), although the difference was not statistically significant (p = 0.0802).

Two case reports [7,8] report a good post-implant outcome without any infection and improved wound healing through the metabolite taurine during post-implant recovery.

PubMed search (<https://pubmed.ncbi.nlm.nih.gov/>): keywords, Tauropace; search date, 06/04/2024.

Ongoing clinical trials

The search on <https://clinicaltrials.gov/> conducted on 24 April 2024 with the keyword 'tauropace' selected the following studies:

- European TauroPace Registry - a University Hospital Schleswig-Holstein Prospective Observational Clinical Study (ETPR) _ NCT04735666. This is an observational, prospective, multicentre study that will evaluate the ability of the TauroPace solution to reduce CIED infections in the first three months after any CIED-related surgical procedure (replacement, revision or de novo implant). The study plans to enrol 2,300 patients and is estimated to end in January 2030.
- Tauropace to Prevent Cardiac Implantable Device Infections in Heart Failure Patients (TRICE-HF) NCT05103267. The primary aim of the multicentre observational study is to evaluate the ability of Tauropace to reduce major CIED infections over 12 months following the CIED procedure in participants with heart failure. The secondary endpoint is to prospectively characterise the performance of Tauropace in participants whose CIED system includes an RV transvenous defibrillation lead. The study will be conducted in up to 10 centres in Italy with up to 500 subjects enrolled. The end of the study is estimated for June 2024.

Guidelines

The latest ESC guidelines [3] recommend a number of pre-, peri- and/or post-procedural actions for the prevention of CIED-related infections, among which are highlighted:

- antibiotic prophylaxis (with flucloxacillin or cefazolin or vancomycin);
- surgical skin preparation (washing the skin with chlorhexidine solution);
- the use of an antibacterial wrap (TYRX) recommended only in high-risk situations.

The guidelines make no reference to the use of disinfectant solutions such as TauroPace to be applied to the CIED prior to implantation.

In a British Heart Rhythm Society (BHRS) document [9] published in January 2024 regarding standards for CIED implantation and follow-up, Tauropace is mentioned as a possible treatment to prevent CIED implant infections in patients at high risk of infection.

Cost-effectiveness analysis

Cost-effectiveness analyses are not available.

MEDLINE search PubMed version conducted on 23 May 2022 with the keyword "(cost[ti] OR economic[ti]) AND (taurolidine OR tauropace)".

HTA Report

Not available

Expected benefits

Reduction of CIED-related infections.

Price and cost of therapy per patient

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)	Annual requirements (No. pieces)
TauroPace (TAUROPHARM GMBH)	200	200	400

Price and cost of therapy per patient with existing treatment alternatives

Product (Manufacturer)	Unit price (euro)	Therapy cost per patient (euro)
------------------------	-------------------	---------------------------------

Price and cost of therapy per patient with existing treatment alternatives

TYRX*	1.114	1.114
-------	-------	-------

* TYRX is not available for purchase to date.

Economic and organisational impact

The incremental cost of post-implant CIED infections is related to the increase of hospital days and ex-implant procedures. According to a retrospective study on a German registry, the incremental costs for post-implantation infection have a range between about 33,000 Euro and about 59,000 Euro [10]. The use of TauroPace could reduce these costs.

Reimbursability information

Main diagnosis (ICD9-CM code)	Intervention (ICD9-CM Code)	DRG code	Tariff (euro)
426 (Conduction disorders)	37.87 (Insertion or replacement of any type of pacemaker)	118 (Cardiac pacemaker replacement)	3.232
		551 (Implantation of permanent cardiac pacemaker with major cardiovascular diagnosis or automatic defibrillator (AICD) or pulse generator)	9.384
		552 (Other permanent cardiac pacemaker implant without major cardiovascular diagnosis)	4.756
		515 (Implantation of cardiac defibrillator without catheterisation)	16.573
		535 (Implantation of cardiac defibrillator with cardiac catheterisation with acute myocardial infarction, heart failure or shock)	24.098
		536 (Implantation of cardiac defibrillator with cardiac catheterisation without acute myocardial infarction, heart failure or shock)	21.634

Assessment of innovativeness (according to Regional Resolution N° 737/2022, link and N° 1244/2022 link)	
Innovative device (Y/N)	N
If yes, please indicate which Criteria 1, 2 and 3 are met:	

Summary data		
Request number	Date of request	Applicant
317	05/04/2024	UOC Cardiology Arezzo- UOSVD Interventional Arrhythmology
Technology		
<p>TauroPace is an extracorporeal disinfectant solution indicated for the eradication of environmental microbial contamination from the surface of CIEDs. Its application is intended for any type of CIED, whether titanium or stainless steel, polyurethane or silicone, including its components (cables, leads, etc.).</p>		
Conclusions		
<p>Implantable cardiac electronic devices (CIEDs) are potentially life-saving treatments for several cardiac conditions, however, they are not without risks. In particular, CIED-related infections are a major issue [1,2]. Despite the prevalence of recommended infection prevention strategies, such as the administration of antibiotics prior to implantation, CIED-related infections continue to increase resulting in increased costs for their treatment [1,3]. Therefore, there is a need to investigate and implement new strategies to reduce the risk of such infections. This includes the TauroPace device, a taurolidine-based disinfectant solution for use in eradicating environmental microbial contamination from the surface of CIEDs prior to implantation in adult patients.</p> <p>Since the first HTA evaluation of TauroPace, a study has become available that compared TauroPace to a solution of H₂O₂ which reported that at 3 months TauroPace reduced CIED infections compared to the control group. The results were significant in that none of the patients assigned to the taurolidine group had an infection during the 3 months following the CIED implant procedure, despite the fact that this group had more risk factors for infection [6].</p> <p>According to the editorial commenting the study [11], the fact that the infection rate with taurolidine treatment was so low is promising. However, more robust, randomised studies are needed to confirm the results of the Borov et al. trial [6]. TauroPace has been mentioned in the British Heart Rhythm Society (BHRS) Guidelines as a possible CIED infection prevention option [9].</p> <p>Overall, the evidence is still limited and there are no direct comparison studies confirming the non-inferiority or superiority of TauroPace over comparators. In spite of these limitations, the OPC expresses a favourable opinion.</p>		
Date of report		
14 May 2024		
Author(s) of the report		
Maria Rita Romeo, Andrea Messori, Sabrina Trippoli		
Referring pharmacist		
Fabiola Del Santo		

BIBLIOGRAPHY

1. Kusumoto FM, Schoenfeld MH, Wilkoff BL, Berul CI, Birgersdotter-Green UM, Carrillo R et al. 2017 HRS expert consensus statement on cardiovascular implantable electronic device lead management and extraction. *Heart Rhythm* 2017;14:e503-51.
2. Dai M, Cai C, Vaibhav V, Sohail MR, Hayes DL, Hodge DO et al. Trends of cardiovascular implantable electronic device infection over 3 decades: a population-based study. *JACC Clin Electrophysiol* 2019;5:1071-80.
3. Blomstrom-Lundqvist C, Ostrowska B. Prevention of cardiac implantable electronic device infections: guidelines and conventional prophylaxis. *Europace*. 2021 May 25;23(Suppl 4):iv11-9. doi: 10.1093/europace/euab071. Epub ahead of print. PMID: 34037227; PMCID: PMC8221047.
4. Shariff N, Akthar T, Razak E, Segerson N, Schwartzman D. Cardiovascular Implantable Electronic Device Infections: Risk Scoring and Role of Antibiotic Envelope in Prevention. *Recent Adv Cardiovasc Drug Discov*. 2015;10(1):70-6. doi: 10.2174/1574890110666151104120207. PMID: 26530457.
5. Clémenty N, Carion PL, Léotoing L, Lamarsalle L, Wilquin-Bequet F, Brown B, Verhees KJP, Fernandes J, Deharo JC. Infections and associated costs following cardiovascular implantable electronic device implantations: a nationwide cohort study. *Europace*. 2018 Dec 1;20(12):1974-1980. doi: 10.1093/europace/eux387. PMID: 29672690.
6. Borov S, Baldauf B, Henke J, Pavaci H, Perani A, Zrenner B, Dietl J, Mehilli J, Lau EW, Vonthein R, Bonnemeier H. Use of a taurolidine containing antimicrobial wash to reduce cardiac implantable electronic device infection. *Europace*. 2023 Oct 5;25(10):euad306. doi: 10.1093/europace/euad306. PMID: 37831737; PMCID: PMC10616572.
7. Giudice M, Catuzzo B, Berlier N, Lau EW, Bonnemeier H, Assadian O, Baldauf B, Borov S, Scacciarella P. Use of Taurolidine in a Patient With a Cardiac Implantable Electronic Device Protrusion. *JACC Case Rep*. 2023 Apr 19;14:101835. doi: 10.1016/j.jaccas.2023.101835. PMID: 37152697; PMCID: PMC10157148.
8. Giaccardi M, Baldauf B, Lau EW, Borov S, Bonnemeier H. Salvage of Cardiac Implantable Electronic Device Pocket Infection with Skin Erosion in Frail 92-Year-Olds. *J Cardiovasc Dev Dis*. 2022 Mar 10;9(3):81. doi: 10.3390/jcdd9030081. PMID: 35323629; PMCID: PMC8955956.
9. British Heart Rhythm Society. Standards for implantation and follow-up of cardiac rhythm management devices in adults. January 2024 revision, <https://bhrc.com/wp-content/uploads/2024/01/BHRS-standards-January-2024-Implantation-and-Follow-Up-of-CRM-Devices-in-Adults.pdf>
10. Ludwig S, Theis C, Brown B, Witthohn A, Lux W, Goette A. Incidence and costs of cardiac device infections: retrospective analysis using German health claims data. *J Comp Eff Res*. 2018 May;7(5):483-492. doi: 10.2217/ce-2017-0080. Epub 2017 Nov 14. PMID: 29132224.
11. Love CJ. Taurolidine: could this be the 'silver bullet' against cardiac implantable electronic device infection? *Europace*. 2023 Nov 2;25(11):euad332. doi: 10.1093/europace/euad332. Erratum in: *Europace*. 2023 Dec 6;25(12): PMID: 37967367; PMCID: PMC10651068.

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Reporting by the Operations Centre, Regional Decree No. 17610 of 7 September 2022.

For further information write to centro.operativo.htart@regione.toscana.it