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Design: A value-based medicine, second-eye model, cost-utility analysis was performed to quantify the comparative effectiveness and cost-effectiveness of therapy with the telescope prosthesis. Participants: Published, evidence-based data from the IMT002 Study Group clinical trial. Ophthalmic utilities were obtained from a validated cohort of >1000 patients with ocular diseases.

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Main outcome measures: Comparative effectiveness was quantified using the (1) quality-adjusted life-year (QALY) gain and (2) percent human value gain (improvement in quality of life). The QALY gain was integrated with incremental costs into the cost-utility ratio (\$/QALY, or US dollars expended per QALY gained).

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4. Hudson HL, Stulting RD, Heier JS, Lane SS, Chang DF, Singerman LJ, Bradford CA, Leonard RE; IMT002 Study Group. Implantable telescope for end-stage age-related macular degeneration: long-term visual acuity and safety outcomes. *Am J Ophthalmol.* 2008 Nov;146(5):664-673. doi: 10.1016/j.ajo.2008.07.003. Epub 2008 Aug 30. PMID: 18760765.
5. Boyer D, Freund KB, Regillo C, Levy MH, Garg S. Long-term (60-month) results for the implantable miniature telescope: efficacy and safety outcomes stratified by age in patients with end-stage age-related macular degeneration. *Clin Ophthalmol.* 2015;9:1099-107.
6. Brown GC, Brown MM, Lieske HB, Lieske PA, Brown KS, Lane SS. Comparative effectiveness and cost-effectiveness of the implantable miniature telescope. *Ophthalmology.* 2011;118:1834-43.
7. Grzybowski A, Wang J, Mao F, Wang D, Wang N. Intraocular vision-improving devices in age-related macular degeneration. *Ann Transl Med.* 2020;8:1549
8. Romano MR, Sorrentino T, Allegrini D, Tripepi D, Stradiotto E, Raimondi R. Small-Incision new generation implantable miniature telescope surgical technique and patient outcomes. *Eur J Ophthalmol.* 2024 Feb 26:11206721241236515. doi: 10.1177/11206721241236515. Epub ahead of print. PMID: 38409807.
9. Savastano A, Ferrara S, Sasso P, Savastano MC, Crincoli E, Caporossi T, De Vico U, Vidal Aroca F, Francione G, Sammarco L, Gambini G, Fedeli C, Di Nardo E, Rizzo S. Smaller-Incision new-generation implantable miniature telescope: Three-months follow-up study. *Eur J Ophthalmol.* 2023 Nov 3:11206721231212545. doi: 10.1177/11206721231212545. Epub ahead of print. PMID: 37920982.
10. Toro MD, Vidal-Aroca F, Montemagni M, Xompero C, Fioretto G, Costagliola C. Three-Month Safety and Efficacy Outcomes for the Smaller-Incision New-Generation Implantable Miniature Telescope (SING IMT™). *J Clin Med.* 2023 Jan 8;12(2):518. doi: 10.3390/jcm12020518. PMID: 36675446; PMCID: PMC9865521.
11. Sasso P, Savastano A, Vidal-Aroca F, Minnella AM, Francione G, Sammarco L, Cima V, Ghirdelli R, Mattei R, Rizzo S. Enhancing the Functional Performance of Patients with Late-Stage Age-Related Macular Degeneration Implanted with a Miniature Telescope using Rehabilitation Training. *Ophthalmol Ther.* 2024 Mar;13(3):697-707. doi: 10.1007/s40123-023-00871-1. Epub 2024 Jan 2. Erratum in: *Ophthalmol Ther.* 2024 Feb 5;: PMID: 38165600; PMCID: PMC10853143.
12. EQ-5D index Calculator, url https://www.economicsnetwork.ac.uk/health/EQ_5D_index_calculator.xls
13. EQ-5D, Versione Italiana, <https://acst-2.org/onewebmedia/italyclin.doc>
14. Balestroni G, Bertolotti G. L'EuroQoL-5D (EQ-5D): uno strumento per la misura della qualità della vita. *Monaldi Arch Chest Dis.* 2012; 78: 155-159

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RAPID HTA REPORT		
No. of request	Date of request	Applicant
315	27-10-2023 (received CO on 15-04-2024)	SOC Ophthalmology at P.O. Palagi (ASL Toscana Centro)
Type of report		
New report		YES
Update of a previous report		
If update, state the reason:		

General technology data					
Trade name					
SING IMT (Smaller Incision New Generation Implantable Miniature Telescope), Model NG SI IMT 3X					
Generic name					
Miniaturized implantable telescope					
Manufacturer name					
Samsara Vision					
Supplier name					
Samsara Vision					
RDM	REF				
2184059	Not reported				
Type	CE mark (date)	Risk class	FDA approval		
Assembled	May 24, 2020	IIB	Yes		
CND					
Q0299 - devices for ophthalmology- Others					
Clinical problem and rationale for the request					
<p>The clinical problem is patients with advanced age-related macular degeneration who are no longer responsive to previous drug treatments. The SING IMT system (second-generation NG SI IMT 3X model) is a device that includes an implantable visual prosthetic telescope and a sterile disposable telescope delivery system for inserting the implant into the eye. The device is designed to improve vision in patients with end-stage age-related macular degeneration and is indicated for bilateral central scotomas.</p>					
Indications for use					
<p>The SING IMT device is indicated for bilateral central scotomas due to end-stage age-related macular degeneration in patients with moderate to profound stable visual impairment.</p> <p>Patients should :</p> <ul style="list-style-type: none"> - Being 55 years of age or older. - Have retinal findings of geographic atrophy or disciform scar with foveal involvement as determined by Angio OCT or fluorescein angiography. - Having signs of cataracts. 					

- Have a distance corrected visual acuity (BCDVA) no better than 20/80 (6/24 metric, 0.25 decimal, 0.60 LogMAR) and no worse than 20/800 (6/240 metric, 0.025 decimal, 1.60 LogMAR) in both eyes.
- Having adequate peripheral vision in the eye not scheduled for surgery.
- Get at least a five-letter improvement on the ETDRS chart in the eye scheduled for surgery when using Samsara Vision's 3X External Telescope Simulator (ETS, supplied separately).
- Have an anterior chamber depth of at least 2.5 mm in the eye scheduled for surgery.
- Be willing to participate in a postoperative training program for the use of the SING IMT.

Patients who have already had cataract surgery are not eligible for SING IMT implantation.

Target patient

Patients with marked central low vision due to end-stage age-related macular degeneration who have not yet had cataract surgery in the eye in which SING IMT is implanted. The SING IMT surgical procedure allows a 2.7-fold magnification of projected images in the part of the retina that is still sensitive, improving visus capabilities.

Main competitors

Currently, these patients have no therapeutic option that can improve central vision impaired by advanced-grade age-related macular degeneration. It should be noted, however, that in a less advanced stage of the disease, some so-called intra-ocular lenses (IOLs) find use, including IOL-AMD system, EyeMax Mono, Scharioth Macula Lens, and Lentis Max, which are described in the review by Borkenstein et al [1].

Technology details

Description

The SING IMTTM system, model NG SI IMT 3X, is an implantable device and delivery system that includes:

- An implantable visual prosthetic telescope referred to below as the "NG SI IMT 3X implant."
- A single-use, sterile telescope delivery system for inserting the implant into the eye, hereafter referred to as the "IMT Delivery System."

The NG SI IMT 3X implant is an implantable prosthetic visual device that, when combined with corneal optical function, is a telephoto system for improving visual acuity in patients with end-stage bilateral age-related macular degeneration.

The NG SI IMT 3X implant includes a glass telescope and haptic mount. The optical component of the telescope, which contains two micro-lenses, is designed to reduce the effective size of the patient's scotoma by magnifying objects in the central visual field and projecting them onto the retina, thus enabling the patient to recognize and identify objects that otherwise could not be seen either from a distance or up close. The optical component (glass telescope) is embedded in a silicon holder, which, after removal of the ocular lens, allows injection of the device into the capsular sac of the eye as part of a conventional ambulatory surgical procedure. The telescope is held in place in the ocular capsular sac by silicone haptic fins.

The telescope delivery system consists of a sterile, single-use cartridge and an injector. The SING IMT cartridge comes pre-loaded with the NG SI IMT 3X implant and injector tip in place ready for loading.

The NG SI IMT 3X implant is surgically implanted in the capsular sac of the eye after lens removal and is held in place by silicone haptic interfaces. Once placed, the NG SI IMT 3X implant renders magnified retinal images of objects in the patient's central visual field.

The NG SI IMT 3X implant is optimized for intermediate distance vision (3 to 10 meters). Corrections for distance and near vision are made with conventional glasses. The device enables viewing of reading materials and other images using natural eye movements rather than head movements. The device provides sufficient image resolution for various activities such as reading, face recognition and TV viewing, with graduated fields of view up to a maximum of 12° (36.0° on the retina). Because there is no relative movement between the eye and the telescope, there are no optical aberrations and a wider field of view is achieved (nominal field of view of 20° projected at about 54° on the retina). Placement of the telescope entirely within the eye eliminates increased motion velocity and vestibular conflict.

The present model is an evolution of the previous one and is characterized by a smaller engraving; it is therefore referred to as the second generation.

Elements of innovation

The first experiences with miniaturized implantable telescope date back many years [2-4]. By querying the PubMed database in combination with the filter on Clinical Trials, the three cited articles [2-4], published between 2006 and 2008, are indeed identified. Other clinical trials, due to obvious keyword insufficiency, are not identified at this stage.

Miniaturized implantable telescopes, developed successively first with a first-generation device (IMT) and then with a second-generation device (which is the subject of this HTA report). The review by Borkenstein et al. [1] is a valuable compendium of evidence on these two devices: the first-generation one was the subject of extensive case histories studied over a long follow-up (more than 200 eyes, studied until more than 5 years [5]), thanks to which the FDA approved the device. The cost-effectiveness was found to be acceptable (\$14,389/QALY), even though the analysis is now many years old [6].

Regarding the second-generation device (referred to as smaller incision new-generation IMT, i.e., SING IMT, the one evaluated here), the evidence is smaller in quantity but more recent [7-10]. The largest case series studied 80 eyes in total [9].

Clinical and economic evidence

Clinical studies

The abstract of the study by Savastano A et al. is quoted in full [9]:

Purpose: To evaluate three months follow-up of SING IMT implant in patients affected by late-stage AMD.

Design: Prospective cohort study.

Subjects: In a total of 80 eyes of 40 patients who underwent the enrollment tests, 11 patients' eyes affected by late-stage AMD matched the inclusion criteria and underwent SING IMT implant from February to June 2022.

Methods: Before surgery, each patient underwent the enrollment examination to verify inclusion and exclusion criteria.

Main outcome measures: BCVA for distance and for near, IOP, ACD and ECD were evaluated at 1 and 3 months follow up. Also quality of life in doing the activities of daily life was evaluated.

Results: BCVA for distance and for near improved from baseline to 3 months follow up (23.91 ± 9.418 ETDRS letters and 59.09 ± 11.58 ETDRS letters respectively ($p < 0.001$). An endothelial cell loss was shown ($p < 0.001$), with a rate of cell density reduction around 8.3% (baseline vs 3 months).

Conclusions: SING IMT could be a valid surgical device to improve patients' sight and quality of life which have been deteriorated by late-stage macular degeneration. Further studies with more patients and longer follow up are needed to confirm our results.

Regarding device-related adverse effects and complications, the most common risks of SING IMT surgery include inflammatory deposits or precipitates on the device and increased intraocular pressure. Significant adverse events include corneal edema, corneal edema that impairs vision, corneal transplantation, and decreased visual acuity. There is a risk that telescope implant surgery may worsen his vision rather than improve it. Individual results may vary.

Ongoing clinical trials

On clinicaltrials.gov (<https://classic.clinicaltrials.gov/ct2/results?cond=&term=sing+telescope&cntry=&state=&city=&dist=>) there are two ongoing studies that are already recruiting patients:

- a) the first is a multicenter trial coordinated by Loma Linda University in California (Concerto trial, NCT05438732). The study plans to enroll 125 patients, and the trial is scheduled to end in June 2024.
- b) The second is a post-marketing observational study (NCT04796545) conducted in Germany. The study plans to enroll 76 patients, and the trial is scheduled to end in February 2027.

Guidelines

Not found.

Cost-effectiveness analysis

The abstract of the U.S. study by Brown et al. [6] is given in full; note that the gain calculated over a 12-year follow-up is 0.76 QALYs.

Objective: To assess the preference-based comparative effectiveness (human value gain) and the cost-utility (cost-effectiveness) of a telescope prosthesis (implantable miniature telescope) for the treatment of end-stage, age-related macular degeneration (AMD).

Design: A value-based medicine, second-eye model, cost-utility analysis was performed to quantify the comparative effectiveness and cost-effectiveness of therapy with the telescope prosthesis.

Participants: Published, evidence-based data from the IMT002 Study Group clinical trial. Ophthalmic utilities were obtained from a validated cohort of >1000 patients with ocular diseases.

Methods: Comparative effectiveness data were converted from visual acuity to utility (value-based) format. The incremental costs (Medicare) of therapy versus no therapy were integrated with the value gain conferred by the telescope prosthesis to assess its average cost-utility. The incremental value gains and incremental costs of therapy referent to (1) a fellow eye cohort and (2) a fellow eye cohort of those who underwent intra-study cataract surgery were integrated in incremental cost-utility analyses. All value outcomes and costs were discounted at a 3% annual rate, as per the Panel on Cost-Effectiveness in Health and Medicine.

Main outcome measures: Comparative effectiveness was quantified using the (1) quality-adjusted life-year (QALY) gain and (2) percent human value gain (improvement in quality of life). The QALY gain was integrated with incremental costs into the cost-utility ratio (\$/QALY, or U.S. dollars expended per QALY gained).

Results: The mean, discounted QALY gain associated with use of the telescope prosthesis over 12 years was 0.7577. When the QALY loss of 0.0004 attributable to the adverse events was factored into the model, the final QALY gain was 0.7573. This resulted in a 12.5% quality of life gain for the average patient during the 12 years of the model. The average cost-utility versus no therapy for use of the telescope prosthesis was \$14389/QALY. The incremental cost-utility referent to control fellow eyes was \$14063/QALY, whereas the incremental cost-utility referent to fellow eyes that underwent intra-study cataract surgery was \$11805/QALY.

Conclusions: Therapy with the telescope prosthesis considerably improves quality of life and at the same time is cost-effective by conventional standards.

HTA Report

Not found

Expected benefits

Improved vision in the target population.

Price and cost of therapy per patient

Product (Manufacturer)	Unit price (euro)	Cost of therapy per patient (euro)	Annual requirement (No. pieces)
SING IMT (Smaller Incision New Generation Implantable Miniature Telescope), Model NG SI IMT 3X (Samsara Vision)	20.000	20,000 plus planting and monitoring expenses	10

Price and cost of therapy per patient with treatment alternatives already in use

Product (Manufacturer)	Unit price (euro)	Cost of therapy per patient (euro)
No therapeutic alternative		

Economic and organizational impact

No specific data are available.

Reimursability information

Principal diagnosis (ICD9-CM code).	Intervention (ICD9-CM Code)	DRG Code	Rate (euro)
362.50 (age-related macular degeneration of the retina - Unspecified)	13.91 (intraocular telescopic lens implantation, telescope implantation)	039 (surgery on the lens, with or without vitrectomy)	039 for 994 euros if DH (refers to cataract surgery) 13.71 for 900 euros if outpatient (refers to cataract surgery)

Evaluation of innovativeness (according to Regional Resolution No. 737/2022, [link](#) and No. 1244/2022 [link](#))

Innovative device (Y/N)	NO
If yes, please indicate which Criteria 1, 2 and 3 are met:	
NOTE: It is emphasized that although a gain of 0.76 QALYs/patient is reported over a 12-year follow-up [6], well above the threshold of 0.15 QALYs posed by the Tuscany Region resolution, these data are derived from simulation analyses (which are frequently the tool for calculating cost-effectiveness) and not from surveys conducted in a "real" setting.	

Summary data

Request number	Date of request	Applicant
315	27-10-2023 (received CO on 15-04-2024)	SOC Ophthalmology at P.O. Palagi (ASL Toscana Centro)
Technology		
SING IMT (Smaller Incision New Generation Implantable Miniature Telescope), Model NG SI IMT 3X, Miniaturized Implantable Telescope.		

Summary data

Conclusions

The device results in a net increase in the cost of patient management, adding to the cost of therapy without replacing a therapy already in use for these cases. At the same time, a net benefit is determined precisely because these patients do not have a therapeutic alternative. The cost-effectiveness is likely to be acceptable even if the analysis simulated the benefit up to a 12-year follow-up, which is far longer than the duration of the trials (3 months in the most recent trial, 5 years in the previous one). In addition, since the analysis was conducted in the United States, the results are not fully transferable to our health care service. Comparative studies are small in number and the duration of their follow-up is still short. The multicenter study, mentioned above, coordinated by Loma Linda University in California, is ongoing.

A detailed proposal considered by the Operations Center is as follows: five treatments at Palagi Hospital managed under the following conditions are initiated on an experimental basis:

- 1) each of these 5 treatments is initially paid for at full price;
- 2) the outcome of these 5 treatments is evaluated, patient by patient, at a follow-up of at least 3 months;
- 3) each patient is assessed as a "favorable outcome" provided that an improvement result of at least 8 letters is achieved (according to the endpoint reported as Visual Acuity Outcomes in the article by Sasso et al [11]); any patient who achieves an improvement of at least 6 letters associated with an increase in utility of at least 0.15 points is also assessed as a favorable outcome (measuring the latter by an after-vs-before comparison using the Italian-language EuroQOL questionnaire [12-14]);
- 4) For patients not classified as "favorable outcome," there is full payback of the cost of the device by the Company.

The Operations Center gives a favorable opinion provided that the points listed above are met. Otherwise, the opinion is understood to be unfavorable.

Date of report writing

May 15, 2024

Author(s) of the report

Andrea Messori

Referring pharmacist

Sara Bellugi

BIBLIOGRAPHY

1. Borkenstein AF, Borkenstein EM, Augustin AJ. Implantable vision-enhancing devices and postoperative rehabilitation in advanced age-related macular degeneration. *Eye (Lond)*. 2023 Mar;37(4):597-606. doi: 10.1038/s41433-022-02179-z. Epub 2022 Jul 22. PMID: 35869389; PMCID: PMC9998400.
2. Colby KA, Chang DF, Stulting RD, Lane SS. Surgical placement of an optical prosthetic device for end-stage macular degeneration: the implantable miniature telescope. *Arch Ophthalmol*. 2007 Aug;125(8):1118-21. doi: 10.1001/archophth.125.8.1118. PMID: 17698761.
3. Hudson HL, Lane SS, Heier JS, Stulting RD, Singerman L, Lichter PR, Sternberg P, Chang DF; IMT-002 Study Group. Implantable miniature telescope for the treatment of visual acuity loss resulting from end-stage age-related macular degeneration: 1-year results. *Ophthalmology*. 2006 Nov;113(11):1987-2001. doi: 10.1016/j.ophtha.2006.07.010. Epub 2006 Sep 20. PMID: 16989902.

4. Hudson HL, Stulting RD, Heier JS, Lane SS, Chang DF, Singerman LJ, Bradford CA, Leonard RE; IMT002 Study Group. Implantable telescope for end-stage age-related macular degeneration: long-term visual acuity and safety outcomes. *Am J Ophthalmol.* 2008 Nov;146(5):664-673. doi: 10.1016/j.ajo.2008.07.003. Epub 2008 Aug 30. PMID: 18760765.
5. Boyer D, Freund KB, Regillo C, Levy MH, Garg S. Long-term (60-month) results for the implantable miniature telescope: efficacy and safety outcomes stratified by age in patients with end-stage age-related macular degeneration. *Clin Ophthalmol.* 2015;9:1099-107.
6. Brown GC, Brown MM, Lieske HB, Lieske PA, Brown KS, Lane SS. Comparative effectiveness and cost-effectiveness of the implantable miniature telescope. *Ophthalmology.* 2011;118:1834-43.
7. Grzybowski A, Wang J, Mao F, Wang D, Wang N. Intraocular vision-improving devices in age-related macular degeneration. *Ann Transl Med.* 2020;8:1549
8. Romano MR, Sorrentino T, Allegrini D, Tripepi D, Stradiotto E, Raimondi R. Small-Incision new generation implantable miniature telescope surgical technique and patient outcomes. *Eur J Ophthalmol.* 2024 Feb 26:11206721241236515. doi: 10.1177/11206721241236515. Epub ahead of print. PMID: 38409807.
9. Savastano A, Ferrara S, Sasso P, Savastano MC, Crincoli E, Caporossi T, De Vico U, Vidal Aroca F, Francione G, Sammarco L, Gambini G, Fedeli C, Di Nardo E, Rizzo S. Smaller-Incision new-generation implantable miniature telescope: Three-months follow-up study. *Eur J Ophthalmol.* 2023 Nov 3:11206721231212545. doi: 10.1177/11206721231212545. Epub ahead of print. PMID: 37920982.
10. Toro MD, Vidal-Aroca F, Montemagni M, Xompero C, Fioretto G, Costagliola C. Three-Month Safety and Efficacy Outcomes for the Smaller-Incision New-Generation Implantable Miniature Telescope (SING IMT™). *J Clin Med.* 2023 Jan 8;12(2):518. doi: 10.3390/jcm12020518. PMID: 36675446; PMCID: PMC9865521.
11. Sasso P, Savastano A, Vidal-Aroca F, Minnella AM, Francione G, Sammarco L, Cima V, Ghirdelli R, Mattei R, Rizzo S. Enhancing the Functional Performance of Patients with Late-Stage Age-Related Macular Degeneration Implanted with a Miniature Telescope using Rehabilitation Training. *Ophthalmol Ther.* 2024 Mar;13(3):697-707. doi: 10.1007/s40123-023-00871-1. Epub 2024 Jan 2. Erratum in: *Ophthalmol Ther.* 2024 Feb 5;: PMID: 38165600; PMCID: PMC10853143.
12. EQ-5D index Calculator, url https://www.economicsnetwork.ac.uk/health/EQ_5D_index_calculator.xls
13. EQ-5D, Italian Version, <https://acst-2.org/onewebmedia/italyclin.doc>
14. Balestroni G, Bertolotti G. The EuroQoL-5D (EQ-5D): a tool for measuring quality of life. *Monaldi Arch Chest Dis.* 2012; 78: 155-159

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Centro operativo

A copy of the document can be downloaded from <http://www.regione.toscana.it/-/prodotti-hta>.

Drafting of the report by the Operations Center, Regional Decree No.17610 of September 7, 2022.

For more information email centro.operativo.htart@regione.toscana.it